



## Stryker 1027 Trauma Stretcher - Fully Refurbished

Stryker's trauma stretcher provides a secure platform for assessing and tending to critically injured and immobile patients. Its innovative design allows for full-length X-ray capability without having to transfer or reposition the patient.

### Features:

- 3" (8 cm) mattress 8" (20 cm) Omni Surface™ casters
- 26" (66 cm) wide patient surface
- Dual cylinder pneumatic backrest
- Dual mounted collection device hooks
- Dual pedestal hydraulics
- Dual side-mounted foot control hydraulics
- Fibresin pneumatic backrest/stationary foot-end litter
- Fold-down siderails
- Four IV receptacles
- Four-wheel, steel-ring brakes with dual end activators
- Full-length X-ray cassette tray
- Integrated storage tray with oxygen bottle holder
- Roller bumpers
- X-ray centering label
- Constant descent hydraulics
- Retractable fifth-wheel steering



### Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.