



Stryker 1550 Electric Stretcher - Fully Refurbished

The Stryker 1550 Electric Stretcher combines the functionality of a hospital bed with the mobility and durability of a stretcher. The versatility of manual to electric operation gives the caregiver the flexibility to handle any situation as patient needs change. Key comfort feature includes electric side-rail controls that allow the patient adjustments. The Stryker 1550 Stretcher has electric knee and back. You can use the select patient controlled adjustments to make it like an electric hospital bed or the manual controls to use it as a standard gurney.

Features:

- Electric controls and manual cranks for head and knee
- Four IV pole receptacles for removable IV poles
- Central locking four-wheel steel ring brakes with dual end activators
- Fifth-wheel steering for enhanced mobility
- Overall Length: 83.5"
- Overall Width: 34"
- Weight Capacity: 500 lbs.
- Height Range High: 35"
- Height Range Low: 22"
- Litter Positioning: Electric
- Brackets: 0-70 deg. Manual
- Backrest: 0-90 deg.
- Knee Gatch: 0-30 deg.
- Trend/reverse Trend: +/- 18 deg.
- Patient Surface: 30" X 75.5"
- Side rails: 13" X 55"
- Caster Diameter: 8"



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.