



## Stryker 721 Transport Stretcher - Fully Refurbished

A variety of available configuration options give this Stryker 721 Big Wheel the versatility you need. The 500 lb weight capacity combined with multiple mobility solutions help accommodate heavy patients and long transports. Positioning options give more flexibility and control than ever to both nurse and patient addressing patient comfort and skin issues associated with extended lengths of stay, helping maintain a positive patient experience.

### Features:

- Pneumatic Fowler Backrest
- Trend and reverse trend
- 75in mattress width (mattress not included)
- Foot pedal position control
- Wheels 20.3(8) Diameter, cm (in)
- Brakes - Central (4)
- Folding siderails
- Steel frame composition
- Steel patient platform
- Big Wheel Steering
- Four IV pole receptacles
- Integrated storage tray with oxygen bottle holder

### Specs:

- Overall length: 83" (210.8 cm)
- Overall width: 30" (76.2 cm)
- Height overall – high position: 36" (91.4 cm)
- Height overall – low position: 24 1/2" (62.2 cm)
- Patient surface: 24" x 76" (61 cm x 193 cm)
- Siderails: 13" x 55" (33 cm x 139.7 cm)
- Caster size: 8" (20.3 cm)
- Backrest: 0° – 90°
- Trendelenburg/Reverse Trendelenburg:  $\pm 18^\circ$
- Maximum patient weight: 500 lb (226.8 kg)



### Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.

