



## Stryker M-Series 104 - Fully Refurbished

The innovative Stryker M-Series SM104 stretchers are ideal for safe transport and effective treatment. Being fully equipped and user-friendly, these stretchers are nicely suited for fast paced and rigorous environments. Completely adaptable, the SM104 is a suitable mobility solution as they reduce the physical strain on clinicians while maintaining patient care and comfort. With the addition of Stryker's BackSmart ergonomic design concepts, caregivers experience less stress-related injury which allows them to improve the patient experience.

### Features:

- Thermoformed ABS base with extra capacity storage tray
- 700 lbs weight capacity
- 26" or 30" (66 cm) wide patient surface with Comfort Mattress
- Less than 21" low height
- Integrated utility tray
- Litter-mounted oxygen bottle holder
- Integrated transfer board
- Power-washable design (SM104 & SM204 only)
- Four-wheel, steel-ring brake system with dual-end activators
- Dual side-mounted foot control hydraulics with uni-lower pedal
- Glideaway™ siderails with integrated push/pull handles
- Four integrated IV receptacles
- Trend./reverse Trend.
- 8" (20 cm) Omni Surface™ casters
- Pneumatic backrest/stationary foot
- Constant descent hydraulics
- Dual-end brake & steer controls
- 5th wheel for additional steering



### Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.