



Stryker Renaissance 1210 Emergency 5th Wheel- Fully Refurbished

The Stryker Renaissance™ Stretcher is a full-featured, durable patient handling system. It is easily maneuverable with standard fifth-wheel steering and pop-up handles. Steel-ring brakes with optional four-sided brake controls ensure a firmly-locked working surface. The Stryker Renaissance 1210 Stretcher can be ordered in various color options and the Big Wheel® option is also available.

Features:

- Height adjustable
- Hand Free Trendelenburg
- 3" (8 cm) mattress
- 8" (20 cm) Omni Surface casters
- 24" (61 cm) wide patient surface
- Fold-to-foot siderails
- Four IV receptacles
- Four-wheel, steel-ring brakes with dual end activators
- Additional 5th Wheel Steering
- Integrated storage tray with oxygen bottle holder
- Locking caster steering
- Roller bumpers
- Stationary foot-end litter
- Push Handles
- IV Pole
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Dual side-mounted foot control hydraulics



Specs:

- Overall Length 84" (213 cm)
- Overall Width 31.5" (80 cm)
- Weight Capacity 500 lb (228 kg)
- Overall Height 22" (56 cm) to 35.5" (90 cm)

Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.