



Stryker SM304 Zoom Series Stretcher - Fully Refurbished

Fully equipped and completely adaptable, the M-Series allows you to choose the mobility solution that best suits your needs. Supremely mobile and user-friendly, these stretchers greatly reduce the physical strain on clinicians without sacrificing patient care and comfort. The Stryker M-Series Stretcher model number SM304 features all standard features plus Zoom Motorized Drive System.

Features:

- Power-Washable design
- Litter-Mounted Oxygen Bottle Holder helps eliminate caregiver bending
- Glideaway™ side rails create a Zero-Transfer Gap
- Integrated transfer board
- Four-wheel, steel-ring brake system with dual-end activators
- Dual side-mounted foot control hydraulics with uni-lower pedal
- Four integrated IV receptacles
- Pneumatic backrest/stationary foot
- Constant descent hydraulics
- Dual-end brake & steer controls

Specifications:

- Overall Width: 37" (94 cm)
- Maximum Patient Weight: 700 lbs (317 kg)
- Overall Size (fully extended length): 85" (215 cm)
- Patient Surface: 30" x 74.25" (76 x 189 cm)
- Max backrest elevation: 0-90 degrees
- Caster Diameter: 8" (20 cm)
- Max knee gatch: 0-40 degrees
- Trend/Reverse Trend ± 15 degrees



Additional Options: Additional IV Poles, Scale, Restraint Straps, Upgraded Mattress Pad

Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rebuild or replace all worn parts or cylinders, and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.