



Tuttnauer 1730 Autoclave - Fully Refurbished

The Valueklave 1730 will satisfy all of your basic sterilization needs at an economical price. The compact size (7" chamber) makes it a perfect primary unit for small offices or a second autoclave for larger practices. The fast cycle time and heated drying delivers a quick turnaround for wrapped or unwrapped instruments. The simple controls and easy closing door make this unit a breeze to operate. High speed and reliability are built into every Valueklave 1730.

Features:

- Fast Cycle Time
- Automatic shut off at the end of both the sterilization and dry cycles
- A long life electro polished chamber and door
- Double safety locking device prevents door from opening while chamber is pressurized
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir
- Dual safety thermostat to protect against overheating



Specs:

- Chamber Dimensions: 7" x 13" (170 mm x 340 mm)
- Chamber Volume: 2 gallons (7.5 L)
- Overall Dimensions: 17.9" (D) x 17.4" (W) x 12.0" (H) (455 mm x 440 mm x 305 mm)
- Tray Dimensions: 11.6" (D) x 4.7" (W) x 0.8" (H) (295 mm x 120 mm x 20 mm)
- Number of Trays: 3
- Standard Unwrapped Cycle Time:
 - Cold: 16 minutes
 - Hot: 11 minutes
- Voltage: 120 V
- Frequency: 50/60 Hz
- Power: 1350 W
- Current: 11 A
- Weight: 50 lbs.

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.