



Tuttnauer 2340M Autoclave - Fully Refurbished

The Tuttnauer manual autoclave is designed to compliment any healthcare facility. The bright new easy to use panel was designed with the operator in mind. Our 2340M autoclave will satisfy all of your sterilization needs without compromising on quality, safety or reliability. This machine utilizes a 316L Stainless Steel chamber and has a dual safety thermostat to protect against overheating. The automatic shut off at the end of the sterilization and dry cycle is an added convenience for you and an additional safety feature for your facility.

Features:

- Automatic shut off at the end of both the sterilization and dry cycles
- A long life electro polished chamber and door
- Double safety locking device prevents door from opening while chamber is pressurized
- Overheat cutoff prevents burnout of heating elements in case of insufficient water in the chamber
- Drain valve is located on the front
- Dual safety thermostat to protect against overheating
- Front panel controls

Specs:

- Chamber Dimensions: 9" x 18.5"
- Volume: 5 gal
- Overall Dimensions D x W x H: 21.5" x 20" x 14.4"
- Standard Cassettes: 2 full/2 half
- Tray Dimensions D x W x H: 6.7" x 16.3" x 0.8"
- No. of Trays: 3
- Standard Unwrapped Cycle Time
 - Cold: 30 mins
 - Hot: 20 mins
- Voltage (V): 120V
- Freq. (Hz.): 50/60
- Power (W): 1400W
- Current (A): 12A



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.