



## Tuttnauer 3870E- Fully Refurbished

The Tuttnauer Model 3870E autoclave is designed for sterilizing laboratory and medical equipment, media and even liquids. The Tuttnauer autoclave provides fast, safe and dependable sterilization without any need for special plumbing or venting. Tuttnauer autoclaves are manufactured with 316L stainless steel sterilization chambers which are highly corrosion-resistant. These autoclaves feature a fast exhaust option which reduces cycle time for solid goods and a slow exhaust option which releases chamber pressure slowly and prevents boiling over of liquids. Cold air pockets are eliminated by continuous pressure purging which ensures uniform sterilization.

### Features:

- Automatic, Digital Control
- Capacity: 84-Liters
- Includes support stand and 2 Stainless Steel Tray Shelves (14 x 26.1-Inch)
- Includes: 60-Minute Timer]
- 3.25-Liter Reservoir & Low Water-Level Alarm
- Temperature Range: 105 to 137°C (221 to 279°F)
- Digital display of temperature & pressure
- Heating time to 273 °F (134 °C): 30-Minutes
- Heat-insulated, double-locking door protects against accidental burns. Door cannot be opened while chamber is under pressure
- Automatic low water level shut-off prevents damage to heating elements when there is insufficient water for the sterilization cycle



### Specs:

- Overall Size: 26"(W) x 34-1/2(L) x 20-3/4(H)
- Chamber: 15" (Diameter) x 30" (L)
- Unit Weight: 220 lbs (99.8 kg)
- Large Tray Size: 14"(W) x 26"(L) x 1"(H)
- Small Tray Size: 11"(W) x 26"(W) 1" (H)

### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.