



Tuttnaur® EZ10 Autoclave - Fully Refurbished

With the simplicity of one touch design all your sterilization and drying needs are fulfilled. The Tuttnauer fully automatic autoclaves fill, sterilize, exhaust and dries at the touch of a button.

Features:

- Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches.
- High efficiency air pump circulates hot air after the sterilization cycle for fast drying.
- Hepa filter 0.2µm air filter (replaceable).
- Pre-loaded cycles for: unwrapped instruments; wrapped/packs; liquid and additional drying cycle.
- All program parameters can be customized and stored.
- Automatic shut off at the end of both the sterilization and dry cycles.
- Audible and visual cycle interruption alert.
- Double safety locking device prevents door from opening while chamber is pressurized.
- Low water sensor prevents activation of sterilization cycle when there is insufficient water in the chamber or reservoir.
- Constant monitoring of temperature and pressure for sterilization.
- Power outage recovery system.
- Door design prevents steam from coming in contact with control panel.



Specifications:

- Chamber Dimensions: 10" x 19"
- Volume: 6 gal / 23 L
- Overall Dimensions: 21.5" x 20" x 14.4"
- Standard Cassette Capacity: 3 Full & 3 Half
- Tray Dimensions: 6.7" x 16.3" x 0.8"
- No. of Trays: 4
- Standard Unwrapped Cycle Time: Cold 30 Min, Hot 20 Min
- Voltage/Frequency: 120 V & 50/60 Hz
- Power/Current: 1400 W / 12 A

Why an Auxo Medical refurbished autoclave?

Auxo autoclaves are refurbished to the highest standards in the industry. We completely disassemble each autoclave to its core and completely rework or replace all parts and return the autoclave to OEM standards functionally and cosmetically. Each autoclave is reassembled, tested and is warranted against defects.