



UMF 4040-650 Power Hi-Lo Exam Table - Fully Refurbished

UMF Standard Series Power Hi-Lo Examination Tables accommodate more patients with a weight capacity of 650 lbs and a lowered height of just 18.5 inches. The Exam Tables also facilitate cleaning with seamless antimicrobial upholstery

Features:

- Power Hi-Lo with Power back adjustment
- Low access ADA compliant 18.5" ($\pm .05$ ")
- 650 lb (294.84 kg) patient weight capacity
- Magnetic hand control with OneTouch Patient Positioning™ two programmable table positioning buttons, plus return-to-home
- Two function foot control with hi-low functions
- Large front storage drawer and reversible rear storage drawer
- Dual electro-mechanical actuator system for efficient speed
- 12 gage steel extension leg-pull out to accommodate for all patient heights
- Removable upholstered top with seamless antimicrobial guard
- Seamless upholstered leg cushion elevates 10° for patient comfort
- Paper roll holder (18" - 21")
- One touch, easy glide stirrups with 4 lateral positions and friction lock length adjustments
- Stainless steel drain pan (1-1/2" deep)
- Hospital grade receptacle 110V with child safety cover (patient left)
- Pelvic tilt (7°)



Specs:

Height: Min. 18.5" ($\pm .5$ ") 46.9cm (± 1.27 cm) at Seat Max. Max: 37" at seat (93.9cm)

Width: 28" (71.12cm)

Length: 57" – 76" with foot section extended (144.7cm – 193cm)

Leg Rest: 19" W x 13" D (48.26cm W x 33cm D)

Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.