



Valley Lab Force Optimum Smoke Evacuator- Fully Refurbished

The OptiMumm Smoke Evacuator is the centerpiece of a complete smoke evacuation system when used with the AccuVac Smoke Evacuation attachment, Valleylab's Universal Mounting System, and either the Force FX Electrosurgical Generator or the Force EZ Electrosurgical Generator.

Features:

- **Quiet:** The ValleyLab OptiMumm System allows you to effectively filter and remove surgical smoke without creating additional noise in the OR.
- **ValleyLab RapidVac System:** Synchronize the ValleyLab OptiMumm Smoke Evacuator with your electrosurgical generator, so the smoke evacuator will work automatically when the generator is activated.
- **Efficient:** The ValleyLab OptiMumm Smoke Evacuator is designed specifically for use with the ValleyLab AccuVac Smoke Evacuation attachment to effectively capture and filter smoke directly at the surgical site.
- **Convenient:** The ValleyLab OptiMumm System is designed with an easy to read front panel and a user-friendly interface. The ULPA and prefilter are easy to access and set up.
- **Compact:** The ValleyLab OptiMumm Smoke Evacuator integrates conveniently into the OR. The system fits directly on most Valleylab mounting carts, saving valuable OR space.
- **User-Friendly:** The filter life indicator is based on airflow, not just time, so you know when to replace the ULPA filter.



Specs:

- Dimensions: 9 ¾" (9" without feet) H x 13 3/8" W x 13 ½" D
- Weight: 35 lbs
- Voltage operating range: 100-120 VAC, 60 Hz; current: 6.0 amps, maximum.

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.