



Zimmer ATS 1500- Fully Refurbished

The Zimmer ATS 1500 is an automatic tourniquet which has been found useful in producing a bloodless operating field in a surgical procedure involving extremities. It is indicated to temporarily occlude blood flow in patient's extremities during a surgical procedure on those extremities. The ATS is found useful in the following surgical procedures: Reduction of Certain Fractures, Kirshner Wire Removal, Nerve Injuries, Tendon Repair, Bone Grafts, Total wrist joint replacements, Knee Joint Replacements, Amputations and Replantations.

Features:

- Line Voltage Range: 105-130 VAC, 50/60 Hz
- Line Current: 485 mA RMS Max., 300 mA RMS nominal
- Input Power: 58 W max, 36 W nominal
- Power Cord Type: SJT, AWG 18, 13 ft. 2 in (4.01m) Power Plug Hospital grade, 3 prong straight blade, 15A
- Line Protection: 0.5A, 250V, 5 x 20 mm SLO-BLO fuse
- Cuff Pressure Range: 50-475 mmHg
- Pressure Accuracy: ± 5 mmHg (50-475 mmHg)
- Pressure Regulation: ± 6 mm of set point
- Maximum Pressure: 475 mmHg (Nominal operation), 600 mmHg (internally limited)
- Time Alarm Set Ranges: 0-240 minutes ; 1 minute increments
- Timer Accuracy: 0.25% ± 1 second
- Battery Type: Rechargeable, sealed lead acid, 12 VDC, 2.3 AH
- Battery Discharge Time: Unit will operate on battery power for 15 minutes (minimum with fully charged battery)
- Battery Recharge Time: 24 hours maximum
- Height: 16.13 in (41 cm)
- Width: 8.75 in (22.3 cm)
- Depth: 7.1 in (17.9 cm) (Including 2 in clamp)
- Weight: 19.8 lbs (9 kg)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.