



Zoll R Series Defibrillator, with SP02 - Fully Refurbished

The R Series ALS continues to meet the needs of ALS providers the world over, and additional features and parameters are being added on a regular basis. Wi-Fi capability is now available. With the addition of a wireless card, R Series can download and transmit code data wirelessly with just the push of a button. Also available is Defibrillator Dashboard, an optional software program that can further reduce your total cost of ownership. Defibrillator Dashboard will automatically notify clinical engineering if a defibrillator is not ready, and even allows for total remote troubleshooting and monitoring of the defibrillator fleet.

Features:

- Real-CPR-Help for ALS
- See Thru CPR filter
- CPR Compression Index
- OneStep system
- OneStep CPR cable
- Smart Prompts
- Code Readiness testing system
- Defib Mentor familiarization aid
- Rectilinear biphasic technology
- 3 lead ECG cable
- SPO2
- Unmatched Clinical Excellence
- Comprehensive Readiness Checks
- Intelligent Tools to Support Users
- SurePower Battery Management
Software for Asset Management
- OneStep simplicity creates new level of rapid, easy deployment
- Dimensions: 8.2" (20.8 cm) (H) x 10.5" (26.7 cm) (W) x 12.5" (31.7 cm) (D).
- Weight: 13.6 lbs (6.2 kg) with OneStep™ Cable and SurePower™ battery pack; 15.2 lbs (6.9 kg) with paddles.



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.