



## Blickman Warming Cabinet - Fully Refurbished Models 7921TG and 7921T

MODELS 7921TG AND 7921TS single door warming cabinets are 74-1/2" tall and ideal for warming blankets or solutions. An open and adjustable shelf design adds to flexibility and expanded use. Easy-to-use controls for setting internal temperature are conveniently located above the door. Each unit comes standard with integrated casters for ease of mobility.

### Features:

- Proportional-Integral-Derivative (PID) digital controllers feature tactile membrane (stainless-steel dome), full numeric keypads, over-temp audible and visual warnings as well as automatic cut-off
- Built-in temperature data recorder with real-time clock
- 160 x 180 pixel graphical display with dual color (red/green)
- Provides "at a glance" trend view, showing a digital representation of recent temperature history
- Self-closing, magnetic-gasket, adjustable tension doors
- Temperature control range allows for setting from 90° to 160°. Solutions only warming cabinets can be easily programmed by user to heat only between 90°F - 110°F
- Self-tune temperature feature maintains set-temperature to within +/- 1°
- Set temperature can be protected two ways: a) Key-Lock b) User defined pass code
- Locking doors standard (except 7927TS/TG)
- No set temperature restriction between compartments
- Calibration due date reminder
- Lost pass code view if user defined pass code is forgotten
- One touch displayed unit conversion from °F to °C and vice versa



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.