



## CMP Blanket & Solution Warming Cabinets - Refurbished

---

CMP has been building high quality warming cabinets for over 65 years and has become the industry standard for hospital and emergency room blanket and solution storage. CMP stainless steel Warming Cabinets can be found in almost every hospital in the U.S. and can be configured for your unique requirements.

**Models:** DSW1AE-18, DSW1AE-18 (74"H), DSW1AE-24, DSW1AE-24 (74"H), DSW2AE-1-18, DSW2AE-1-24, DSW2AE-2-18, DSW2AE-2-24, DSW2AE-2-24, DSW3AE-18 (18"W), DSW3AE-24

### Features:

- Stainless Steel Construction
- Electric Heated Compartments
- Over Temperature Alarm
- Separately Controlled Chambers
- Temperature Adjusting Lockout
- Optional "Quick-Temp" USB temperature monitoring system
- Many Standard Sizes Available
- Microprocessor Control System
- Customizable Temperature Settings/Limits
- Wire Shelves/Baskets to maximize airflow
- Removable Pilaster Strips for easy cleaning
- Optional Wi-Fi and Ethernet connectivity
- Optional 3/4" thick Thermopane Glass Doors
- Digital Display
- UL/CUL Listed
- $\pm 2^{\circ}\text{F}$  Accuracy
- Uniform Heating
- Optional Casters



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.



## CMP Warming Cabinet Specifications Per Model



Model	Outside Dimensions	Compartment Dimensions	Volume	Weight	Compartments	Shelves
DSW1AE-18	30" W x 45" H x 21" D (762 x 1143 x 533 mm)	26.56" W x 25.88" H x 17.50" D (675 x 657 x 445 mm)	7.0 cu. ft.	180 lbs (82 kg)	1	1 fixed, 1 adjustable
DSW1AE-18 (74"H)	30" W x 74" H x 21" D (762 x 1880 x 533 mm)	26.56" W x 54.87" H x 17.50" D (675 x 1394 x 445 mm)	14.8 cu. ft.	390 lbs (177 kg)	1	1 fixed, 4 adjustable
DSW1AE-24	30" W x 45" H x 27" D (762 x 1143 x 686 mm)	26.56" W x 25.88" H x 23.50" D (675 x 657 x 597 mm)	9.4 cu. ft.	200 lbs (91 kg)	1	1 fixed, 1 adjustable
DSW1AE-24 (74"H)	30" W x 74" H x 27" D (762 x 1880 x 686 mm)	26.56" W x 54.87" H x 23.50" D (675 x 1394 x 597 mm)	19.8 cu. ft.	410 lbs (186 kg)	1	1 fixed, 4 adjustable
DSW2AE-1-18	30" W x 74" H x 21" D (762 x 1880 x 533 mm)	26.56" W x 21.75" H x 17.50" D (Top) (675 x 552 x 445 mm) 26.56" W x 30" H x 17.50" D (Bottom) (675 x 762 x 445 mm)	13.9 cu. ft.	400 lbs (181 kg)	1 heated, 1 unheated	2 fixed, 3 adjustable
DSW2AE-1-24	30" W x 74" H x 27" D (762 x 1880 x 686 mm)	26.56" W x 21.75" H x 23.50" D (Top) (675 x 552 x 597 mm) 26.56" W x 30" H x 23.50" D (Bottom) (675 x 762 x 597 mm)	18.7 cu. ft.	420 lbs (191 kg)	1 heated, 1 unheated	2 fixed, 3 adjustable
DSW2AE-2-18	30" W x 74" H x 21" D (762 x 1880 x 533 mm)	26.56" W x 21.75" H x 17.50" D (Top) (675 x 552 x 445 mm) 26.56" W x 30" H x 17.50" D (Bottom) (675 x 762 x 445 mm)	13.9 cu. ft.	400 lbs (181 kg)	2	2 fixed, 3 adjustable
DSW2AE-2-24	30" W x 74" H x 27" D (762 x 1880 x 686 mm)	26.56" W x 21.75" H x 23.50" D (Top) (675 x 552 x 597 mm) 26.56" W x 30" H x 23.50" D (Bottom) (675 x 762 x 597 mm)	18.7 cu. ft.	420 lbs (191 kg)	2	2 fixed, 3 adjustable
DSW3AE-18	30" W x 26" H x 21" D (762 x 660 x 533 mm)	26.56" W x 14.25" H x 17.50" D (675 x 362 x 445 mm)	3.8 cu. ft.	160 lbs (73 kg)	1	1 fixed
DSW3AE-18 (18"W)	18" W x 26" H x 21" D (457 x 660 x 533 mm)	14.56" W x 14.25" H x 17.50" D (370 x 362 x 445 mm)	2.1 cu. ft.	140 lbs (64 kg)	1	1 fixed
DSW3AE-24	30" W x 26" H x 27" D (762 x 660 x 686 mm)	26.56" W x 14.25" H x 23.50" D (675 x 362 x 597 mm)	5.2 cu. ft.	180 lbs (82 kg)	1	1 fixed