



auxoMEDICAL

Family Practice Exam Table

For comfort, economy, overall functional and storage, there's nothing better than Clinton's Family Practice Table. This versatile table has storage access from either side. Plus, its lower treatment table height could eliminate the need for a step stool, and makes it ideal for family practice, clinic or wherever a cost effective treatment/exam table is needed.

Features:

- Premium, vinyl upholstery over 2" firm, foam padding (5 cm)
- Full 27" width for patient comfort (68.58 cm)
- Air-spring adjustable backrest
- Coated, two position adjustable stirrups
- Pull-out leg rest on steel track
- Optional pull-out basin
- Easy care all laminate base
- Color coordinated doors and top available
- Chip resistant edges
- 1 large side storage compartment with shelf
- Storage compartment is accessible from both sides (no right or left handed tables)
- 2 drawers with all steel slides
- Built with plywood and laminate panels that are California CARB II compliant
- Top sections are easily replaceable in the field
- Adjustable feet for leveling
- Paper dispenser included
- 400 lbs. load capacity under normal use (181.43 kg)



Options:

Pull-Out Pan, Paper Cutter

Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.