



Hausted Unicare III 800 Series Stretchers - Fully Refurbished

The 800 series stretcher has a thermoformed base cover made of durable polymer for easy cleaning. The design offers you a stretcher that is lightweight, durable and corrosion resistant.

Features:

- Low Height: (59.7cm)
- High Height: (81.3 cm)
- Patient Surface
 - Regular Width: (64.1 x 193.0 cm)
 - Extra Wide: (76.2 x 193.0 cm)
- Overall Dimensions
 - Regular width: (73.7 x 201.9 cm)
 - Extra Wide: (83.8 x 201.9 cm)
- Side Rails
 - Retracto Rails: (149.9 x 39.3 cm)
 - Vertical Rails: (147.3 x 36.8 cm)
- Casters: (20.3 cm diameter)
- Brake and steer pedals: 4 corners
- IV and accessory wells: 4 corners
- X-Ray capability: Optional
(standard on 829)
- Fowler Options on Pneumatic, Crank and Manual
- Knee Gatch: Optional
- Trendelenburg/Reverse Trendelenburg: 12"
- Easy-lift II pneumatic assisted fowler backrest, stainless steel, operable from both sides of the stretcher.



Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.