



## Hologic Flouroskan Premier Encore Mini C-Arm - Fully Refurbished

Premier Encore builds on a strong foundation of delivered ease of use and time saving functionality. With features like 6/4 Dual Mode and High Speed Mode, you'll discover the perfect balance of precision and productivity — from imaging larger areas like knees and shoulders to the tightest magnification on fingers and toes.

### Features:

- Digital acquisition and a 45-micron focal spot produce high-definition images with superior clarity
- SureLock, a patented safety feature, locks the C-arm securely into position
- The iso-centered C-arm rotates a full 360° and rolls smoothly into position around the surgical subject
- 6/4 Dual Mode lets you switch from a 6-inch (15 cm) to a 4-inch (10 cm) field of view
- 15 inch (38 cm) high intensity monitors swivel 60° for easy viewing at any angle and images are sharp and clearly defined even in the brightest surgical environment
- Touch-sensitive right-hand display allows direct selection of images, operating commands and Mosaic Mode displays
- Moves and maneuvers easily on 6-inch wheels and safety brake locks system firmly in place
- Premier Encore offers more storage options: 2000 or 4000 image hard drive storage, 3.5-inch floppy diskette or 3-inch and 5-inch CD-R/CD-RW or DICOM transfer
- Triple function foot switch includes x-ray on/off, image save and image print



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.