



Medfusion 3500 V. 4 Syringe Infusion Pump - Fully Refurbished

The unique, easy-to-use, and easy-to-learn design of the Medfusion 3500 Syringe Pump with PharmGuard Medication Safety Software gives you the ability to customize the platform to your institution's practices. In other words, it adapts to the way you work, not the other way around.

Features:

- FlowSentry Enhanced pressure trending for earlier detection of occlusion.
- Extensive delivery mode combinations including ml/ hr., body weight, mass, volume over time, custom dilution and intermittent, loading dose, bolus dose, standby, volume limit and KVO rate.
- Syringe sizes: Accepts seven different manufactures' syringe types and all manufactures' syringe administration sets.
- +/- 2% instrument accuracy: Automatic syringe size detection.
- Ergonomic, easy-to-use, horizontal design protects entire syringe barrel and allows single-handed loading.
- MRI indicated: For use outside the 150 Gauss line and with shielded magnets of eld strength of 1.5 Tesla or less.
- Mounts easily to an IV pole, an infant's isolette or radiant warmer. Pole clamp allows for rotation to accommodate various pole and vast orientations.
- Bio-Med mode to set PM dates, view service history, and for access to calibration and diagnostic testing.
- Internal battery capacity of 10 hours at 5.0 mL/hr. with 60 ml syringe
- Wide-ranging low rates and the capacity to infuse at low rates with larger syringes
- Dimensions: 10.5" wide x 5.1" high x 5.5" deep (27 x 13 x 14 cm)
- Weight (pump only): 3.8lbs. (1.7kg)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.