



Medfusion 3500 V6 Syringe Pump - Fully Refurbished

The Medfusion 3500 V6 Syringe Pump with PharmGuard Medication Safety Software is an intuitive system designed to enhance safety. The pump has unlimited configurable clinical alerts supporting best practice at the bedside.

Features:

- Wide screen allows up to 30 characters for each drug name to help reduce the need for abbreviations
- 240 x 64 pixel liquid crystal display. An always on LED driven backlight enhances viewing in low light areas
- Weight, volume and occlusion programmable limits
- Loading and bolus dose therapies
- Alarm volume, times and type
- FlowSentry rapid occlusion detection continuously monitors pressure trending
- Horizontal syringe plunger for easy loading, unloading and barrel protection
- Battery capacity of 10 hours fully charged (infusing 5.0 ml/hr with 60 ml syringe, 25° C)
- Electromagnetic Compatibility:
Tested to EMI/EMC testing per IEC 60601-1 clause 21.6.
- Dosing Modes are user customizable
- Dimensions: 270 w x 130 h x 140 d mm
- Weight: 1.7 kg / 3.8lbs (Pump only)
- Accuracy: +/- 2%
- Delivery Rate Range: 0.01 - 1130 ml/hr
- Minimum Infusion Increment: 0.01 ml/hr



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.