



auxoMEDICAL

Midmark 604 Exam Table - Fully Refurbished

Midmark's 604 Manual Examination Table gives you the best economy exam room option available with choice features. And this cost-effective table does not skimp on patient comfort. With the 604, the physician can easily position patients so that they are comfortable and relaxed throughout the exam.

Features:

- Exam Assistant Fold-down Drawer System provides an area to store basic supplies and can be used as an extra work surface when needed
- Flex drawer design allows changing to either right- or left-hand access all without the use of tools
- Premium Comfort System is designed to enhance patient comfort for a chair-like feel in an exam table
- Large, rubberized footstep with bright yellow edging designed to assist the visually impaired
- Steel-frame construction is designed to support patients up to 500 pounds
- Height with Upholstery: Premium/Ultra-Premium Top: 33.38" (84.8 cm)
- Length with Footrest Extended: 75.5"
- Patient Weight Capacity: 500 lb (226 kg)
- Electrical Receptacle: (115 VAC, 5 AMPS Max.)
- Duplex, hospital grade, electrical receptacle on patient's left side
- Paper Roll Storage: 4 Rolls - 21" x 3.5" (53.3 cm x 8.9 cm)
- Paper Roll Holder: Holds 21" x 3.5" Paper Roll (53.3 cm x 8.9 cm)
- Stirrups are standard



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.