



Mindray DPM 4 Monitor — Fully Refurbished

The DPM 4 monitor is a compact, mobile anesthesia companion made for the versatile clinic. Continuous patient monitoring on a high-resolution color display allows for uninterrupted anesthesia care, no matter how many places it's needed.

High-efficiency is of paramount importance to anesthesia clinicians and the DPM 4 was created for that purpose. Clear, precise data transmission and easy integration with current anesthesia equipment means more patient time and less mechanical configuration.

Features:

- 96 hours of graphic and list trends
- 40 seconds of full disclosure waveform recall
- 3 trace recorder for printing real-time and historical data
- Compact ash back-up
- Nurse call interface with connection to hospital call systems
- 8.4" TFT display
- Lithium ion battery
- 3/5-lead ECG
- Dual-channel invasive blood pressure
- Non-invasive blood pressure
- Nellcor®
- OxiMax®
- SpO2
- Pulse rate respiration
- Dual-channel temperature
- Lithium ion battery
- Integrated recorder



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.