



auxoMEDICAL

NUVO VisitOR LED Surgical Lights - New

The NUVO Surgical Light provides superior performance for all aspects of surgery. The light head provides excellent shadow reduction, a high-intensity light pattern and optimal color rendition. Its unique, patented interface allows the user to fully control the light's intensity, pattern, position and orientation from the light head. Appropriate for the vast majority of general and specialized surgical procedures including Micro Surgery, Orthopedics, OB/GYN, Ophthalmology, and Cosmetic Surgery.

Features:

- Light Intensity (LUX): 120,000
- Color Temperature Adjustment: 4300°K
- Pattern Intensity: 5 Settings, 50 -100 %
- Pattern Size (Inches): 9-13
- Red Rendering Index: 95
- Controls in the Sterile Field: Yes
- LED Projected Life (hrs): 50,000

Available:

- Single Light Head
- Single Light Head + Monitor
- Single Light Head + Camera
- Dual Light Head
- Dual Light Head + Monitor
- Dual Light Head + Camera
- Triple Light Head
- Triple Light Head + Monitor
- Triple Light Head + Camera



Why purchase equipment from Auxo Medical?

At Auxo Medical, we don't take any decision lightly. Whether we focus on quality of new or refurbished medical equipment, repair savings, or improved operational cost-efficiency—our eye is always on the smallest details that can affect the level of care provided and your budget. Our consultants pride themselves in helping you plan and develop the most appropriate facility for your changing needs.