

Ritter Midmark 354 Light - Fully Refurbished

Whether you practice in an office environment, specialized procedures area or a hospital emergency room, your crucial work requires dependable, quality lighting. The kind of lighting you get from Midmark. The Ritter 354 Procedure Spotlight has everything healthcare professionals need - but without the heat. It's an ideal task light for all medical environments where illumination is required for external examinations and procedures. And with its pure, crisp and defined light pattern, it can be adjusted to suit individual needs.

Features:

- Illumination at 24" (60.9 cm)
- Focal length: 7,000 fc (75,300 lux)
- Illumination at 36" (91.4 cm)
- Focal length: 4,000 fc (43,000 lux)
- Maximum pattern diameter at 36" (91.4 cm)
- Focal length: 10" (25.4 cm)
- Minimum pattern diameter at 36" (91.4 cm)
- Focal length: 3" (7.6 cm)
- Color temperature: 4,200 K Diameter of lighthead: 5" (12.7 cm)
- Bulb: One 21 VAC, 150-watt halogen lamp
- Lighthead pivot: 180° Vertical axis pivot: 570°
- Counterbalanced for precise positioning
- Available in single and dual mount and mobile version
- Includes an easy to adjust dimmer lever that controls the amount of light, and another level which adjusts the spot size from 3" to 10".
- Built-in fan on the back of the light head pulls airflow through the light to keep cool.
- Fully autoclavable handle



Why purchase Auxo Medical Refurbished Lights?

Auxo Lights, halogen and LED, are refurbished to the highest standards in the industry. We completely disassemble each light and arm assembly, and completely rebuild or replace all worn parts and return each assembly to OEM standards. Each Light is reassembled, tested and is warranted against defects.