



Ritter Midmark 355 Procedure Light - Fully Refurbished

Whether you practice in an office environment, specialized procedures area or a hospital emergency room, your crucial work requires dependable, quality lighting. The kind of lighting you get from Midmark. The Ritter 355 is the perfect balance between quality and affordability. It's counterbalanced to provide precise positioning without drifting. A 4,200 K color temperature, large pattern size with an even distribution of light, and unsurpassed shadow control provides the most accurate color rendition for improved diagnosis.

Features:

- Illumination at 24" (60.9 cm)
- Focal length: 5,000 fc (53,800 lux)
- Illumination at 36" (91.4 cm)
Focal length: 3,400 fc (36,600 lux)
- Pattern diameter: 8"
- Depth of field: 35"
- Color temperature: 4,200 K
- Diameter of lighthouse: 17"
- Bulb: One 100-watt halogen lamp
- Requires less than 2 ½ pounds of force to move the light
- Fully autoclavable handle
- Truer color temperature with less heat
- Illuminated 8-inch diameter pattern size
- Available in single or dual ceiling mounts as well as five-caster mobile base models



Why purchase Auxo Medical Refurbished Lights?

Auxo Lights, halogen and LED, are refurbished to the highest standards in the industry. We completely disassemble each light and arm assembly, and completely rebuild or replace all worn parts and return each assembly to OEM standards. Each Light is reassembled, tested and is warranted against defects.