



auxoMEDICAL

Skytron Integrity 215 Sterilizer - Fully Refurbished

Integrity 215 exceeds all standards for capacity—efficient enough for immediate use, productive enough for high-demand environments, and compact enough for most replacement applications. Building on Skytron's history of quality and innovation, the 215 is part of a family of sterilizers designed to meet the wide-ranging needs of today's healthcare environments. Whether you're outfitting a new facility, renovating an older one, or replacing outdated technology, Integrity sterilizers offer industry-changing efficiency, capacity, ease of use, and use of space.

Features:

- Spacious 21 1/2" x 21 1/2" x 38" sterilization chamber
- Processes up to six 25-pound trays (150 pounds per cycle)
- Integrated water recirculation system
- 5.7" LCD full-color touchscreen panel
- Countdown timer and progression bar visible from more than 30'
- 12 standard and 12 fully customizable cycles
- Electric-powered sliding door
- Recessed foot-pedal and touchscreen door control
- Stainless steel and copper piping
- Dual-walled, fully jacketed stainless steel pressure vessel
- Programmed automatic cycle start
- Energy saving daily startup/shutdown timer
- Electronic recording of up to 9,999 cycles
- Integrated ink-to-paper printer for hard copy cycle reports
- Password protected access to supervisor and service modes



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.