



auxoMEDICAL

Skytron Integrity 270 Sterilizer - NEW

Integrity 270 is the workhorse of the Integrity line, exceeding all standards for capacity and processing up to sixteen 25-pound trays per cycle. Building on Skytron's history of quality and innovation, the Integrity 270 is part of a family of sterilizers designed to meet the wide-ranging needs of today's healthcare environments. Whether you're outfitting a new facility, renovating an older one, or replacing outdated technology, Integrity sterilizers offer industry-changing efficiency, capacity, ease of use, and use of space.

Features:

- Processes up to sixteen 25-pound trays (400 pounds per cycle)
- Integrated water recirculation system
- 5.7" LCD full-color touchscreen panel
- Countdown timer and progression bar visible from more than 30'
- 12 standard and 12 fully customizable cycles
- Electric-powered door locking/unlocking mechanism
- Stainless steel and copper piping
- Dual-walled, fully jacketed stainless steel pressure vessel
- Electronic recording of up to 9,999 cycles
- Integrated ink-to-paper printer for hard copy cycle reports
- Programmed automatic cycle start
- Energy saving daily startup/shutdown timer
- Spacious 26 1/2" x 37 1/2" x 53" sterilization chamber



Why purchase equipment from Auxo Medical?

At Auxo Medical, we don't take any decision lightly. Whether we focus on quality of new or refurbished medical equipment, repair savings, or improved operational cost-efficiency—our eye is always on the smallest details that can affect the level of care provided and your budget. Our consultants pride themselves in helping you plan and develop the most appropriate facility for your changing needs.