



auxoMEDICAL

Skytron J2 Series Warming Cabinets - Fully Refurbished

Skytron's J2 Series Warming Cabinets are designed to maintain clinically ideal temperatures for solution bottles, solution bags and blankets. Six models are available to provide optimum clinical storage requirements.

Counter Top Models: SS2201-J2, SS2201-J2G

Free Standing Models: SS2205-J2, SS2205-J2G, SS2207-J2, SS2207J2G

Features:

- Double wall stainless steel construction with number 4 finish
- Field reversible left or right hinged doors
- Optional double paned, tempered glass doors
- Heating compartments fully insulated with circulating fan
- Two or three adjustable stainless shelves
- Digital temperature controls for individual compartments
- Recirculating fan in each compartment
- Adjustable temperature range
- Set point temperature +/- 5° F
- Independent power switch / compartment
- Audible and visual over-temperature indicator
- Automatic shut off (temperature exceeding 180° F)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.