

Steris AMSCO QDJ06 Warming Cabinets -Refurbished

The Steris AMSCO Dual Compartment Warming Cabinet is designed to address high-volume patient warming needs. This model provides up to 15.8 cu. ft. capacity and offers the versatility of two independently controlled warming compartments. Designed to be more than a blanket warmer, the Dual Compartment Warming Cabinet can store and warm sterile intravenous fluids (IV fluids), surgical irrigation fluids, and linens and blankets.

Features:

- LED readout for instant temperature confirmation
- Superior durable materials to withstand years of use
- Adjustable shelving for flexible organization and storage
- Reliable heating with uniform and efficient heat transfer
- Visual access to inventory with dual-pane tempered glass or traditional solid steel doors
- Glass or stainless steel door and adjustable shelving
- Temperature lock-out for consistency and control
- Celsius and Fahrenheit selection for convenience and accuracy
- Outside dimensions (H x W x D) at 18" Depth: 74 3/4 x 30 x 20 1/2"
- Outside dimensions (H x W x D) at 24" Depth: 74 3/4 x 30 x 26 1/2"
- Total Capacity: at 18" depth = 11.6 cu. ft.
- Total Capacity: at 24" depth = 15.8 cu. ft.
- Temperature Range: 90°F (32°C) to 160°F (71°C)
- Door Option: Steel or dual-pane tempered glass
- Free-standing or recessed installation
- Optional electronic numeric keypad door lock and electronic data reporting with USB port



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.