

System Two LED Series Lights - New

The System Two LED Surgery Lights represent the pinnacle of the System Two series. Utilizing the latest LED technology, these lights are the brightest, longest-lasting light source in our history. Each light head offers 120,000 lux output with 50,000 hour-rated LED life. The 4,300° K light color offers a beautifully white light field that provides excellent visualization. The five-stage dimming control and 9"– 13" light field are both adjusted entirely with the sterilizable handle, putting all control in the sterile field. Available in several different mounting configurations, along with video capability, the System Two LED offers a tremendous value in O.R. lighting.

Features:

- 120,000 lux output at 1 meter
- 4,300° Kelvin color temperature
- 9"– 13" adjustable pattern
- Five-stage dimming
- 50,000 hour-rated LED life
- Ability to control all functions in the sterile field
- Provides options for High Definition, fiber optics cabling, video source switching and compatibility with PACS
- Can be ordered with or without monitor
- Ceiling mount version as solo arm or shared with lights
- Sterile handle for adjustments within the sterile field
- 1080P high-definition video camera
- HD-SDI video output
- 20x optical zoom/12x digital zoom
- Auto and manual focus



Why purchase equipment from Auxo Medical?

At Auxo Medical, we don't take any decision lightly. Whether we focus on quality of new or refurbished medical equipment, repair savings, or improved operational cost-efficiency—our eye is always on the smallest details that can affect the level of care provided and your budget. Our consultants pride themselves in helping you plan and develop the most appropriate facility for your changing needs.