



Durasonic 6.6 Gal Digital Ultrasonic Cleaner- Fully Refurbished

Each model of the Durasonic Digital Ultrasonic Cleaners is constructed using durable industrial style 40kHz transducers. These provide increased cleaning power along with built in sweep frequency to ensure uniform cleaning activity throughout the bath. The five larger sizes have built in drains. Each model has a digital timer and mechanical thermostat, with a power switch/fuse in the rear of the machine.

Features:

- Industrial Grade Model
- Sweep and Degas Function
- Temperature Controlled Heat for best cleaning performance
- Operator Temperature Controls (knob)
- Digital Ultrasonic Timer
- Comes with Free Basket and Cover
- Very Reliable Industrial-Grade Transducer
- Tank Material: Stainless Steel
- Includes faucet-type drain
- 2-year warranty
- Heat Light Indicator
- Timer Control Buttons; set timer only when ultrasonics is off.



Specs:

- Inside Dimensions (LxWxH) : 18.5"x11.5"x8"
- Includes Heater
- Sonic Power Watt Avg: 1000
- Includes Basket
- Weight: 60 lbs.
- Power Supply: AC 100 ~ 120v, 50 / 60Hz
- Fuse Size Amp: 15

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.