



Gendron 1190-231 Bariatric Stretcher- Fully Refurbished

Gendron develops, designs and manufactures bariatric patient care products for acute care, long-term care, home care and rehabilitation environments. The Gendron 1190-231 has a patient weight capacity of up to 1000 pounds and a litter top usable patient surface that measures 37" x 80".

Features:

- Rugged reinforced frame and litter top
- Three position vertical side rails
- Foot pedal control hydraulic fowler back rest
- Four wheel central locking
- Full perimeter non-marring bumper
- Four IV receptacles
- 4" high density mattress pad
- Overall bumper to bumper width 40" for door way access
- Litter top to floor height measures 29 3/8"
- Beachwood powder coat finish
- Patient weight capacity up to 1000 pounds
- Warranty: 1 year parts only warranty



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.