



## Hausted Eye Surgery Chair- Fully Refurbished

The Hausted Surgi Chair is a multi-position chair and stretcher intended for treatment, transport, or recovery of a patient. With the uniquely designed articulating head section, this unit can be used for specialty procedures that require precision head positioning. Models in this series provide specialized surgical functions of head positioning for ophthalmology, plastic surgery, or ENT procedures.

### Features:

- Low Profile Height Adjustable Articulating Head Piece with Dual Operation Control Knobs
- Motorized Controls for Height Adjustment: 22.4" to 32", Back Rest: 0°- 80°, Leg Section
- Trendelenburg: 0°- 6°
- Patient/Staff Hand Held Control
- Emergency Backrest Release Lever to Override Power Backrest
- 24 VDC Battery Back-Up Allows Power Functions When Not Plug Into An Outlet
- Patient Surface: 22" x 75"
- Weight Capacity: 500 lbs
- Customizable colors
- Included Mattress
- Sectional 3" Thick Mattress Pad
- Warranty: 1 year
- Powered controls to easily maneuver head and leg sections for patient comfort.
- Quick positioning with Trendelenburg controls and emergency release turning critical movements.
- With the hand held control, patients can position themselves during recovery, which minimizes staff involvement and can help shorten recovery time.
- Easy accessibility getting on and off chair with a low height of 19".
- Height Range: High 80 cm
- Height Range: Low 60 cm



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.