



Lifepak 12 with CO2, Aux and Pacing- Fully Refurbished

The LIFEPAK® 12 Defibrillator/Monitor is a multi-parameter device that combines semi-automated and manual defibrillation with capnography, external pacing, 12-lead electrocardiography and other monitoring functions. Its platform is upgradeable, allowing customers to add new parameters as their needs change and as new features are introduced.

Features:

- Available with ADAPTIV™ biphasic technology
- Compact and portable
- Compatible with QUIK-COMBO® pacing/defibrillation/ECG electrodes
- With Co2, Aux and Pacing
- Easy, simple, and affordable upgrades
- Masimo SET® pulse oximetry offers accurate and stable oxygen saturation monitoring
- Oscillometric noninvasive blood pressure (NIBP) monitoring
- End-tidal CO2 (EtCO2) monitor for both intubated and nonintubated patients, with patented Microstream®, capnography technology and FilterLine® accessories that operate smoothly even in high humidity
- Low Battery Indication and Message
- 31.7 cm (12.5 in) H x 38.9 cm (15.3 in) W x 21.7 cm (8.5 in) D
- The device is capable of transferring data records by PC Card internal modem, external EIA/TIA modem, cellular modem, or serial connection.
- Heart Rate Display
- Voice Prompts
- Alarms: Quick Set, VF/VT Alarm, Apnea Alarm, Heart Rate Alarm



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.