



Lifepak 20 With Pacing & SPO2- Fully Refurbished

The LIFEPAK® 20 defibrillator/monitor is the ideal crashcart defibrillator putting early, effective defibrillation into the hands of BLS responders while transitioning easily to full therapy and monitoring capabilities for the ALS-trained code team. The 20 is highly intuitive, making it easy for infrequent AED-trained responders to quickly understand and use.

Features:

- Simple yet sophisticated for responders of many skill levels
- Highly intuitive for quick, effective AED use
- Easily converts to manual mode
- ADAPTIV™ biphasic technology adjusts shock duration and voltage
- Broad dosage capability up to 360J, when needed
- Compact and lightweight for easy carrying and versatile use
- Superb color graphics allow easy viewing from many angles
- With Pacing and SpO2



Specs:

- Waveform: Biphasic truncated exponential
- Energy: ± 1 joule or 10% of setting, whichever is greater, into 50 ohms ± 2 joule or 15% of setting, whichever is greater
- Charge Time: 200J < 4 seconds & 360J < 7 seconds with fully charged battery, 360J < 10 seconds while not in low battery operations
- Sweep Speed: 25 mm/sec for ECG
- Frequency Response: 0.67 to 32 Hz
- Heart Rate: 20 - 350 BPM
- Battery: NiMH
- Humidity: 5 to 95%, noncondensing
- Vibration: MIL-STD-810E Method 514.4, cat1

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.