



MR-Conditional Transport Stretcher- Fully Refurbished

MRI Stretchers are extremely beneficial to any MRI suite that requires transportation of prone patients. Due to the inherent power of the magnet in MRI machines, only non-magnetic equipment can be used within the MRI room. MRI stretchers are an important addition for MR facilities that requires patient transport in and around the MRI suite. The MR-Conditional Transport Stretcher is designed for general patient transport, the non-magnetic transport stretcher features all-stainless-steel construction for MRI use.

Features:

- 2" litter pad
- Swing-down side rails
- Full perimeter non-marring bumper
- Manual fowler backrest for patient positioning
- Four individual all-surface locking caster wheels
- Fully welded frame and litter top
- Shipped set up no assembly needed
- Four I.V. receptacles includes one I.V. Pole
- MR-conditional to 3-Tesla
- Load capacity: 350lbs
- Patient surface: 76" L x 23" W
- Warranty: 1 Year



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.