



Sharpertek Ultrasonic Cleaner 17 Gallon- Fully Refurbished

SharperTek® professional ultrasonic cleaning systems meet the demands of industry and commerce. Their attention to detail in the design and manufacturing process and using only the highest quality components results in a professional ultrasonic cleaning system that you, and your business, can count on.

Features:

- Industrial grade model
- Heating element for best cleaning results
- Heating element can stay on all day even if ultrasonic power is off
- Stamina XP Series®
- Digital Ultrasonic Timer
- Industrial-Grade Transducer for best reliability
- Drain Faucet for spill-proof draining
- Very reliable, powerful, industrial-quality transducers, 40 KHz
- Temperature-controlled heating



Specifications:

- Ultrasonic Frequency: 40,000 Hz
- Adjustable feet provide stabilization and prevent rocking
- Tank Material: Stainless Steel
- Very reliable, powerful, industrial-quality transducers, 40 KHz
- 14-Gauge Stainless Steel tank and body
- 17 gallon tank
- Tank dimensions: 22" × 12" × 15" (L×W×H)
- Unit dimensions: 35" × 33" × 35" (L×W×H)
- Temperature-controlled heating
- Power Supply AC 100 ~ 120V, 50 / 60Hz
- Power 2400 W (1200 Ultrasonic and 1200 Heating)
- Weight: 180 Lb
- Warranty: 2 Years

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.