



Steelco US 100L Ultrasonic Cleaner- Fully Refurbished

The Steelco US 100L is a single tank ultrasonic bath utilized for the treatment of reusable medical devices and rigid endoscopes. It combines a stainless steel cleaning tank with industrial style transducers and a powerful ultrasonic generator to provide the strongest table top cleaning available. It combines the speed and effectiveness of ultrasonic cleaning power with the convenience of plug-in-anywhere operation.

Features:

- Stainless steel basin and cabinet stand
- External dimensions 35.5 x 27.5 x 35.5 (WxDxH) (tray lift closed)
- Automated instrument tray lift for hands-free operation / immersion 3
- Level -loading of trays, usable space per level 25" x 13" x 4" (WxDxH)
- Hands-free operation by foot-pedal for opening / closing
- 5 customizable cleaning programs
- LED display indicating temperature and remaining cycle time
- Automatic filling, dosing and drain
- 6 connections for flushing of cannulated instruments
- Stand cabinet with storage area for detergent
- 1 dosing pump with level controller for liquid detergent
- Stainless steel stand with pull-out drawer for detergents
- Built-in ultrasonic generator with adjustable power
- Capacity for 3 DIN instrument trays
- 12 luer lock connections
- 2 removable instrument supports
- USB port
- Automatic minimum level control of chemical with alarm
- Internal rinsing of cannulated instruments
- Built in detergent storage compartment
- Holds up to 19" x 10" x 2" baskets



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.