



Steris Caviwave Single Chamber- Fully Refurbished

The Caviwave Ultrasonic Cleaning Systems are freestanding units that automatically clean, flush, lubricate, and dry instruments (depending on unit configuration) prior to final wash and disinfection cycles*. Provides powerful micromechanical scrubbing frequency of 132kHz which more than triples the commonly used 40kHz.

Features:

- Cleaning frequency of 132kHz that more than triples the commonly used 40kHz
- Toe touch cycle start
- Automatic chemistry injection
- Automatically cleans, flushes, lubricates, and dries depending on unit configuration*
- Flush ports allow tank solution reach into internal passages
- Can process up to 35lbs. of instruments

Sizes Available:

- 15-gal (57 L) tank holds 22 lb (10 kg) of instruments. tray size 20 x 13-1/2 x 3-1/2"
- 20-gal (76 L) tank holds 35 lb (16 kg) of instruments. tank sizes 24 x 13-1/2 x 3-1/2"

Specs:

- 208 V, 50/60 Hz, 3-Phase (5 Wire)
- Drain
 - 1-1/4" Slip CPVC (trap and waste lines are not provided by STERIS).
- Hot Water
 - 1/2" FPT CPVC (38 psi maximum); 110-120°F (43-49°C).
- 39-1/2 x 42-1/4±1 x 30"(1003x1073±25x762mm)
- Service Access Clearance: Minimum 24" (610 mm) all sides
- Warranty: 1 year warranty man. defect, 90 day PM parts



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.