



## Winco Stretchair- Fully Refurbished

The Winco Stretchair is the solution to safer patient transfer; this chair encourages a no-lift culture in most healthcare environments. This versatile chair can go from upright to infinite recline to lateral transfer positioning by using the attached hand-control pendant. Individually adjustable armrests and leg rests/footrests give more comfort for the patient.

### Features:

- Powered adjustment of height and recline position operated with a four-button hand controller.
- Independent leg rest operation with flip-up and adjustable length footrests.
- Seat height adjusts from 22" – 30" (55.8 – 76.2 cm).
- Lay-flat stretcher position adjusts from 30" – 38" (76.2 – 96.5 cm).
- 1000 lb. patient weight capacities.
- Lap/safety belt, integrated oxygen cylinder holder and IV pole with dual position mounts
- Large rear-mounted storage bin.
- Quick connect high-capacity battery.
- On-board battery recharging system with battery charge indicator.
- Coiled power cord easily and quickly stores out of the way.
- Centrally controlled multi-surface 6" (15.25 cm) casters with four-wheel lock and two-wheel steering.
- Full width push handle in seating position for easy maneuvering of chair.
- Powder-coated frame and metal components for long, durable life.
- Available in a variety of designer upholstery colors



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.