



# auxoMEDICAL

## ConMed Sabre Genesis - Fully Refurbished

The Sabre Genesis™ provides state-of-the-art technology in a compact and affordable unit perfect for use in office-based electrosurgical procedures. The Genesis™ incorporates ConMed's Dynamic Response technology, which reacts to different tissue impedances in less than 1 millisecond ensuring a consistent delivery of power and clinical effect while minimizing surrounding tissue damage. The Sabre Genesis™ is versatile enough to handle a variety of office-based procedures with modes that include: pure cut, blend cut, pulse cut, coag, and bipolar.

### Features:

- Dynamic Response Technology: Controls current and voltage at 450,000 times per second for optimal clinical effects
- Continuous ARM Safety Monitoring: Advance warning when dispersive electrode contact is compromised
- Patient safety features reduce blood loss and the risk of accidental burns
- Integrated interface for activation with a smoke evacuator
- Easy set-up
- Combination handswitch/footswitch jack for multiple accessory connections
- Front panel volume control
- Option to select automatic programming which restores the Sabre Genesis™ to last settings used
- Pulse Cut Mode provides precise energy bursts for critical dissections
- Choice of modes and surgical effects ideal for in-office procedures
- Weight: 14.65 lb. (6.65 kg)
- Dimensions: 6" (H) x 13" (W) x 18" (D) including handle



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.