

TABLE OF CONTENTS

Anesthesia/Equipment Booms 2

 Medispace Pendant System Single Retractable 3

 Valia X Equipment Arm with M20 Pendant System Spring Arm System
 3

 Mediboom Pendant System Tandem Mounted Fixture 3

 Amico iCE LP Electric Series - Single Mount Dual Arm 1000mm 4

 Amico Tandem Light and Anesthesia Universal Mount 5

 Nuvo Navigator Pendant System Tandem Mounted Fixture 6



ANESTHESIA/EQUIPMENT BOOMS

Auxo Medical provides anesthesia/equipment booms that are flexible. They conform to various configurations depending on the space requirements, function and purpose. In fact, anesthesia/equipment booms reduce clutter caused by the cords they eliminate the risk of tripping hazards. Ultimately, this prevents hazards, thereby improving staff and patient safety. Anesthesia/Equipment Booms are also a great option for surgical lighting systems, making them an efficient and ideal solution for any facility.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

[REFURBISHED](#)




MEDISPACE PENDANT SYSTEM SINGLE RETRACTABLE

The height-adjustable MediSpace, a ceiling-mounted spring arm system with integrated media connections, reliably and cost-effectively supplies floor-based machines such as anesthesia machines, heart-lung machines or endoscopic equipment trolleys with medical gases, electricity and data. Thanks to the integrated compensating spring, the MediSpace can be precisely positioned in space with one hand and easily brought to an ergonomic working height. In this way, the connection between the medical devices and the MediSpace can be easily established in the lower position. Thanks to the generous travel, the system can be mounted close to the ceiling and, thanks to the space-saving parking position, gives the OR team more freedom of movement compared to conventional anesthesia ceiling supply units.

[Read More](#)

SKU: AM-MEDISPACE

Price: Request Quote



VALIA X EQUIPMENT ARM WITH M20 PENDANT SYSTEM SPRING ARM SYSTEM

VALiA X is a slim support arm system for loads of up to 110 kg for the height-adjustable version with spring arm and up to 140 kg for the non-height-adjustable extension arm version. The combination of the slim design with the high load capacity of the spring arm makes VALiA X the ideal support arm system for large screens in diagnostic rooms. Thanks to its slim design, VALiA X, combined with a supply console, also offers a space-saving solution for intensive care units.

[Read More](#)

SKU: AM-VALIA-X-M20

Price: Request Quote



MEDIBOOM PENDANT SYSTEM TANDEM MOUNTED FIXTURE

The Mediboom has the most versatile array of features in the industry. The combination of duplex outlets and provision for gases combined with the ease with which it can be moved, the Mediboom Equipment Management System is capable of supporting your full range of equipment in today's operating room. The use of tapered roller bearings and a healthy throat size on the suspension arms allow the Mediboom to accommodate a full complement of services without compromising mobility. Service capacity and smooth, easy-to-move arms with rock solid brakes have become trade marks of the Mediboom Equipment Management System.

[Read More](#)

SKU: AM-MB-TAND

Price: Request Quote

AMICO ICE LP ELECTRIC SERIES - SINGLE MOUNT DUAL ARM 1000MM

Each Unit Consists of:

- (1) 18" x 9" C30 Console
- (2) DISS Oxygen Outlet and Hose Assembly for US
- (1) DISS Medical Air Outlet and Hose Assembly for US
- (1) DISS N2O Outlet and Hose Assembly for US
- (1) DISS WAGD Outlet and Hose Assembly for US
- (2) DISS VAC outlet and hose assembly for US
- (4) 20 A Hospital Grade Emergency Duplex Receptacle Assembly (2) Blank Provision
- (2) Blank Provision for Data
- (1) 28x23 LT Shelf with side rail

[Read More](#)

SKU: LP-L-11N-N-E

Price: Request Quote



AMICO TANDEM LIGHT AND ANESTHESIA UNIVERSAL MOUNT

Each Unit Consists of:

- Rods & Fasteners for light Ice 30 Light Fixture
- (1) Electric Brake Power Control Unit (1 Power Supply) (1) Gas Riser Bracket (15 Riser)
- (1) Riser Kit for AIR
- (1) Riser Kit for N2O
- (2) Riser Kit for OXY
- (2) Riser Kit for VAC
- (1) Riser Kit for WAGD

[Read More](#)

SKU: LH-MNT-UAP-DBL-BOM

Price: Request Quote



NUVO NAVIGATOR PENDANT SYSTEM TANDEM MOUNTED FIXTURE



The Navigator™ has the most versatile array of features in the industry. The combination of duplex outlets and provision for gases combined with the ease with which it can be moved, the NUVO Navigator Equipment Management System is capable of supporting your full range of equipment in today's operating room. The use of tapered roller bearings and a healthy throat size on the suspension arms allow the Navigator to accommodate a full complement of services without compromising mobility. Service capacity and smooth, easy-to-move arms with rock solid brakes have become trade marks of the NUVO Navigator Equipment Management System.

Three Sizes of Service Head and Utilities -- Compact (model 413), Intermediate (model 600) and Large (model 1025)

[Read More](#)

SKU: AM-NV-NAV-TAND

Price: Request Quote

INDEX

A

Amico iCE LP Electric Series - Single Mount Dual Arm 1000mm 3

Amico Tandem Light and Anesthesia Universal Mount 4

M

Mediboom Pendant System Tandem Mounted Fixture 3

Medispace Pendant System Single Retractable 2

N

Nuvo Navigator Pendant System Tandem Mounted Fixture 5

V

Valia X Equipment Arm with M20 Pendant System Spring Arm System 2