

TABLE OF CONTENTS

Autoclaves	3
SteriData Logger Paperless Sterilizer Cycle Printer For Statim	4
SteriData Logger Sterilizer Printer, Autoclave Printer, Paperless - Ritter / Midmark M9 / M11 Gen1 Series	4
Tuttnauer T-TOP 10B Class B Pre/Post Vac Autoclave	5
Tuttnauer T-TOP 11B SFPP Autoclave	5
Alternative Data Logger to Midmark 9A682001 for M9/M11-040 through - 042	6
Alternative Data Logger to Midmark 9A682001 for M9/M11-020 through - 022	7
Alternative Data Logger to Midmark 9A682001 for M9/M11-001 through - 019	7
Enbio PRO Class B Autoclave	8
Enbio S Class B Autoclave	8
Tuttnauer T-TOP 11S SFPP Autoclave	9
Tuttnauer 3870M Autoclave (New)	9
Scican Statim 5000 G4 Cassette Autoclave (New)	10
Midmark Ritter M9 UltraClave (Autoclave), M9-022	10
Midmark Ritter M11 UltraClave Veterinary Autoclave, M11-033	11
Midmark M9D Manual Door Dental Autoclave	11
Tuttnauer 2540EK Kwiklave Electronic Autoclave (Refurbished)	12
Cassette for Scican Statim G4 5000 Autoclave (Refurbished)	13
Cassette for Scican Statim G4 5000 Autoclave (New)	13
Steri-Dent 200 Dry Heat Sterilizer	14
Tuttnauer T-TOP 10S Gravity Autoclave	15
Tuttnauer LABSCI 85LV Vertical Autoclave	15
Tuttnauer LABSCI 60LV Vertical Autoclave	16
Tuttnauer LABSCI15L+ Benchtop Autoclave	17
Tuttnauer LABSCI15L Benchtop Autoclave	17
Tuttnauer LABSCI11L Benchtop Autoclave	18
Tuttnauer DS1000 Steam Distiller 1 Gallon	18
Tuttnauer DS1G Steam Distiller 4L - 1L per Hr	19
Tuttnauer 3870HSG-WS Pre/Post Vacuum 85L / 22.5 Gallon Autoclave	19
SteriData Logger Sterilizer Printer, Autoclave Printer, Paperless - Ritter /	

Midmark M9 / M11 NEW Series	20
Tuttnauer 3870E Autoclave	20
Midmark Ritter M11 Autoclave, M11-001	20
Tuttnauer 3870EA Autoclave Gen 2 (New)	21
3870EA Autoclave Gen 2 (Refurbished)	21
T-Edge 10 or 11 Barcode Printer	22
Tuttnauer T-Edge 11S Gravity Autoclave	22
Tuttnauer T-Edge 11B Pre Post Vacuum Autoclave	23
Midmark Ritter M11D Dental Manual Autoclave	24
Tuttnauer EZ 11 PLUS No Printer	24
Tuttnauer EZ10k 220V Autoclave	25
Ritter M11 Veterinary Autoclave	26
Tuttnauer Elara 11 Autoclave	26
Ritter M11 Autoclave, M11-001 - RENTAL	26
Midmark M9/M11 Printer and cable	27
Ritter Midmark M9, M11, Old Style, Connection Cable Included (Retired)	27
Statim 2000 G4 Cassette Sterilizer	28
Tuttnauer 2540M Autoclave	28
Tuttnauer 2540MK	29
Cox 6000 Rapid Dry Heat	29
Dri-Clave 3000 Dry Heat Sterilizer 3 Tray	30
Statim 5000 Sterilizer	30
Tuttnauer 3870M Autoclave (Refurbished)	31
Scican Statim 5000 G4 Cassette Autoclave (Refurbished)	32
Dri-Clave Dry Heat Sterilizer - 2 Tray	32
Amsco Eagle 10 Plus	33
Pelton & Crane Validator 10	33
Midmark Ritter M11 Ultraclave (Autoclave), M11-022	33
Tuttnauer 2340M Manual Autoclave	34
Tuttnauer Valueclave 1730 Autoclave	35
Pelton & Crane Magnaclave	35
Tuttnauer EZ10	35
Ritter Midmark M9, M9-001	36
Tuttnauer 3870EA Autoclave	36
Statim 2000 Sterilizer	37
Ritter M11 Autoclave, M11-001	37
Tuttnauer 2340EA Autoclave	38
Tuttnauer EZ9 Plus	38



AUTOCLAVES

For every need and budget we provide efficient autoclaves for medical facilities and laboratories. Our inventory consist of reliable brands, such as [Tuttnauer](#), [Ritter](#), and more.

Our safe, durable, FDA-approved machines come in a variety of sizes and instrument capacities. Auxo Medical's inventory of autoclaves is reliable and compact enough to fit into whatever space your facility offers.

We offer both small- and large-format autoclaves. Ultimately, we understand the importance of maintaining sterilized instruments for any medical facility.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

[REFURBISHED](#)



STERIDATA LOGGER PAPERLESS STERILIZER CYCLE PRINTER FOR STATIM

SteriData™ Logger for SciCan Statim 2000 and 5000 cassette sterilizers*
(must have must have serial port in back and current software)

The paperless alternative to a paper Printer

This SteriData Logger unit plugs in just like the paper printer and sits on top of your Sterilizer. Now all your sterilizer cycle data is saved to the SD card (included) automatically.

When the cycle data is required, simply take the SD card to your computer where all the cycle files are arranged by year and month, date and time.

The individual sterilizer cycle files are recorded as text files and easily displayed and/or printed out from your computer as required. Optionally, you can include your practice name and the sterilizer serial number which will be added to every cycle file.

[Read More](#)

SKU: AM-S1219DL

Price: \$1,075.00



STERIDATA LOGGER STERILIZER PRINTER, AUTOCLAVE PRINTER, PAPERLESS - RITTER / MIDMARK M9 / M11 GEN1 SERIES

SteriData™ Logger for Midmark Ritter M9 and M11 Ultraclaves, M3 Ultrafast sterilizer, and SciCan Statim 2000 and 5000 cassette sterilizers*

The paperless alternative to a paper Printer

This SteriData Logger unit plugs in just like the paper printer and sits on top of your Sterilizer. Now all your sterilizer cycle data is saved to the SD card (included) automatically.

When the cycle data is required, simply take the SD card to your computer where all the cycle files are arranged by year and month, date and time.

The individual sterilizer cycle files are recorded as text files and easily

displayed and/or printed out from your computer as required. Optionally, you can include your practice name and the sterilizer serial number which will be added to every cycle file.

[Read More](#)

SKU: AM-M1219DL-G1

Price: \$1,095.00



TUTTNAUER T-TOP 10B CLASS B PRE/POST VAC AUTOCLAVE

The **Tuttnauer T-TOP 10B Class B Pre/Post Vac Autoclave** is a compact, highly efficient sterilization solution designed for ease of operation and advanced connectivity. Featuring a 10" W x 18" D chamber with a 21L capacity, it accommodates up to five trays and offers sterilization times ranging from 24 to 73 minutes depending on the cycle, with drying included. Its user-friendly interface, equipped with a multicolor LCD touchscreen, simplifies setup and daily operation, while programmable cycle scheduling saves time. The autoclave's advanced technology includes internal memory capable of storing up to 500 cycles, remote access via Wi-Fi and T-Connect cloud software for cycle monitoring, notifications, and paperless documentation—all supporting streamlined workflow and enhanced safety. Designed for reliable performance, it features steam flush pressure pulse (SFPP) technology without liquids, and offers options for software updates to ensure ongoing efficiency. Its sturdy design and manual water filling make it suitable for a variety of healthcare and sterilization centers seeking dependable, connected sterilization.

[Read More](#)

SKU: AM-T-TOP10B

Price: Request Quote




TUTTNAUER T-TOP 11B SFPP AUTOCLAVE

The **Tuttnauer T-TOP 11B SFPP Autoclave** offers advanced, reliable sterilization with user-friendly features designed to enhance efficiency in healthcare settings. This Class B pre/post vacuum autoclave utilizes steam flush pressure pulse (SFPP) technology, ensuring thorough sterilization without liquids, with programs optimized for unwrapped, wrapped, and delicate instruments—ranging from 38 to 80 minutes at various temperatures. Its 27-liter chamber (11" W x 18" D) accommodates multiple trays, and it is easily operated via a multilingual LCD touchscreen featuring intuitive icons. The autoclave is equipped with internal memory capable of storing up to 500 cycles, programmable start times, and advanced connectivity options including two USB ports, Ethernet, Wi-Fi, and T-Connect cloud-based IoT software—allowing remote access to cycle history, notifications, and real-time monitoring. Designed for modern sterilization centers, it features a paperless documentation system, advanced programming options, and a specially designed user interface to streamline daily workflow, making it a robust, efficient solution for ensuring sterilization compliance and safety.

[Read More](#)

SKU: AM-T-TOP11B

Price: Request Quote




ALTERNATIVE DATA LOGGER TO MIDMARK 9A682001 FOR M9/M11-040 THROUGH -042

The **Autoclave Data Logger M42** is an alternative to the Midmark 9A682001 for third-generation Midmark M9 and M11 sterilizers with LCD displays and gray-white covers, specifically models ending with -040 through -042. Measuring approximately 2-1/2" x 4" x 1", this wireless device includes a removable 16GB micro SD card, SD card adapter, stand, cable, user/setup manual, and a drill template for installation under the printer cover plate. It features a secure Wi-Fi connection, FCC certification, and an LCD display. The M42 automatically logs sterilization cycle data, maintains accurate time with daylight savings adjustments, and uploads logs seamlessly to a secure offsite Dropbox cloud backup. With a five-year battery life (CR1220) and unique device naming, it provides a reliable, offsite data management solution tailored for third-generation Midmark autoclaves.

[Read More](#)

SKU: 9A682001-M42

Price: \$795.00




ALTERNATIVE DATA LOGGER TO MIDMARK 9A682001 FOR M9/M11-020 THROUGH -022

The **Autoclave Data Logger M22** is a reliable alternative to the Midmark 9A682001 data logger for Midmark M9 and M11 sterilizers numbered -020 through -022. Sized at approximately 2-1/2" x 4" x 1", this wireless device includes a removable 16GB micro SD card, SD card adapter, stand, cable, user/setup manual, and a drill template for easy installation under the printer cover plate. Featuring a secure Wi-Fi connection, LCD display, and FCC certification, the M22 automatically logs sterilization cycle data, synchronizes time with daylight savings adjustments, and uploads data seamlessly to your offsite Dropbox cloud backup. With a five-year battery life (CR1220) and unique naming for each autoclave, it provides a simplified, secure solution for offsite data management.

[Read More](#)

SKU: 9A682001-M22

Price: \$795.00



ALTERNATIVE DATA LOGGER TO MIDMARK 9A682001 FOR M9/M11-001 THROUGH -019

The **Autoclave Data Logger M19** serves as an effective alternative to the Midmark 9A682001 data logger for older-style Midmark M9 and M11 sterilizers ending with -001 through -019. Designed for compatibility with Generation 1 Midmark autoclaves featuring a red display, this wireless data recorder measures approximately 2-1/2" x 4" x 1" and includes a removable 16GB micro SD card, SD card adapter, stand, cable, user/setup manual, and a drill template for installation under the printer cover plate. Equipped with a secure, Wi-Fi-enabled LCD display, the M19 automatically logs, stores, and uploads sterilization cycle data to Offsite Dropbox cloud backup, supports automatic time and daylight savings adjustments, and features a 5-year battery life (CR1220). FCC certified and uniquely named for each autoclave, this device offers reliable, offsite data management for your sterilization processes. Please contact us prior to ordering if you have any questions.

[Read More](#)

SKU: 9A682001-M19

Price: \$795.00



ENBIO PRO CLASS B AUTOCLAVE

The **Enbio PRO Autoclave** is a state-of-the-art, FDA-approved Class B pre- and post-vacuum sterilizer designed for medical and dental practices seeking rapid, reliable sterilization in a compact footprint. Weighing only 45 lbs and up to three times smaller and lighter than conventional sterilizers, it delivers cycle times—including drying—in as little as 32 minutes, making it ideal for busy office environments. Capable of sterilizing up to 16 instrument packs and suitable for hollow and complex instruments, it features quick cycles starting at just 9 minutes for unwrapped items at 134°C and comprehensive pre-post vacuum sterilization. Easy to install—simply plug in—and comes complete with user manuals, validation reports, and Enbio Data Viewer software on a USB drive. The waste bottle is sold separately. This innovative sterilizer offers unmatched efficiency, portability, and performance in a sleek, space-saving design.

[Read More](#)

SKU: AM-ENBIO-PRO

Price: Request Quote



ENBIO S CLASS B AUTOCLAVE

The **Enbio S Autoclave** is a compact, lightweight, and highly efficient Class B pre- and post-vacuum sterilizer designed for medical and dental clinics. With cycle times—including drying—of just 32 minutes, it offers rapid, reliable sterilization while occupying up to six times less space and weighing three times less than traditional autoclaves. Its versatile design accommodates hollow instruments with quick start cycles, beginning at just 7 minutes for unwrapped items at 134°C. Easy to operate—simply plug in—and complete with user manuals, validation reports, and Enbio Data Viewer software on a USB stick, the Enbio S delivers FDA-approved performance in a small footprint, making it ideal for office and clinical environments. Note that waste bottles are sold separately.

[Read More](#)

SKU: AM-ENBIO-S

Price: Request Quote



TUTTNAUER T-TOP 11S SFPP AUTOCLAVE

The **Tuttnauer T-TOP 11S SFPP Autoclave**—a compact, high-performance sterilizer designed to simplify and accelerate your daily workflow. Its user-friendly LCD color touchscreen and intuitive icons make programming and monitoring effortless, while a multistage steam flush (SFPP) ensures fast, liquid-free sterilization. With internal memory for up to 500 cycles, remote access to cycle history via Wi-Fi and T-Connect software, and cloud-based paperless documentation, you can manage, customize, and optimize operations from anywhere. The device supports time-saving scheduling, delivers precise cycle options for wrapped and unwrapped loads, delicate instruments, and small quick-turnaround runs, and provides secure cloud connectivity for alerts and notifications. Built for reliability, the 11L chamber (27L total volume) features 5 trays and a compact footprint (23" D x 19.5" W x 17.5" H) with easy water filling and dual USB ports, making it an ideal sterilization centerpiece for busy offices or clinics seeking streamlined, connected, and compliant sterilization.

[Read More](#)

SKU: AM-T-TOP11

Price: Request Quote



TUTTNAUER 3870M AUTOCLAVE (NEW)

The **Tuttnauer 3870M Manual Autoclave** offers reliable, efficient sterilization with user-friendly features like automatic shut-off at the end of both sterilization and drying cycles. Crafted with a durable, electro-polished chamber and door for long-lasting performance, it includes a double safety locking device to prevent accidental opening while pressurized. The front drain valve allows for quick water reservoir draining, and dual safety thermostats provide added protection against overheating. With a spacious 15" x 30" chamber (85L capacity), it accommodates up to 15 full standard cassettes and includes trays measuring 26.6" x 11.3" and 26.6" x 13.8". Designed for efficiency and safety, this autoclave operates on 230V, 50/60Hz power, weighing approximately 209 pounds, making it a dependable choice for medical and laboratory sterilization needs.

[Read More](#)

SKU: AM-3870M-N

Price: Request Quote



SCICAN STATIM 5000 G4 CASSETTE AUTOCLAVE (NEW)

The SciCan Statim 5000 G4 Cassette Autoclave delivers a fully automatic and exceptionally fast sterilization solution for today's demanding practices. This brand-new unit utilizes four single-touch programs to optimize instrument cleaning from start to finish. Boasting a larger cassette and an additional sterilization cycle, the Statim 5000 G4 efficiently processes even larger instruments with remarkable speed and reliability, completing a wrapped cycle in just 17.5 minutes and an unwrapped cycle in a mere 9 minutes. The modernized design features a user-friendly 3.5" color touchscreen with intuitive controls, and a built-in data logger that eliminates the need for traditional thermal paper and cumbersome external devices by saving cycle data to a USB drive or the cloud. Equipped with features like ultra-fast cycles, Dri-Tec Drying Technology, and aseptic transfer capabilities, the Statim 5000 G4 ensures exceptional instrument turnover, prolongs instrument life, and reduces instrument investment, making it an invaluable addition to any practice.

[Read More](#)

SKU: AM-STATIM-5000-G4-N

Price: Request Quote



MIDMARK RITTER M9 ULTRACLAVE (AUTOCLAVE), M9-022

The **Midmark Ritter M9 UltraClave** is a user-friendly automatic sterilizer designed for simple and reliable operation. Its automated features streamline the sterilization process: simply load the instruments, select the desired cycle, and press start. The M9 automatically fills with the correct water level and, upon cycle completion, automatically opens the door for efficient steam dissipation and faster drying times. Built-in safety features, including overheat protection, a pressure relief valve, and an error light, ensure safe and reliable performance. Self-program controls allow customization of cycle parameters to accommodate specific instrument processing needs, while pre-programmed cycles offer convenience for routine sterilization. The green LCD display provides clear visibility of cycle indicators. The M9 is commonly used in various healthcare settings such as dental clinics, medical offices, veterinary practices, and other outpatient facilities where efficient and reliable sterilization of instruments is critical. Its ease of use and compact design make it well-suited for smaller practices and labs. Industries beyond healthcare, such as tattoo parlors and beauty salons, may also utilize the M9 for sterilizing tools and equipment.

[Download Brochure](#)

[Read More](#)

SKU: AM-M9-022

Price: Request Quote



MIDMARK RITTER M11 ULTRACLAVE VETERINARY AUTOCLAVE, M11-033

The **Midmark Ritter M11 UltraClave Veterinary Autoclave**, including the refurbished model with an 11" x 18" chamber, is renowned for its remarkable sterilizing capacity, making it an ideal choice for veterinary practices that handle larger packs and instruments. This countertop sterilizer boasts an exclusive automatic-opening door, which not only enhances ease of use but also facilitates a quick drying time, ensuring efficiency in every cycle. Its Steam-Flush Pressure Pulse Air Removal System guarantees proper steam penetration, thereby optimizing sterilization effectiveness. Users can take advantage of the Self Program Controls to tailor cycle parameters for specific instrument processing needs, while the Pre-Programmed Controls offer added convenience for routine tasks. The green LCD display provides a clear view of operational indicators, enhancing user interaction. With external dimensions of 17.8" wide by 23.8" deep by 17.8" high, the unit accommodates two trays measuring 9" x 15" and two trays measuring 6.6" x 15". Weighing 90 pounds, the M11 delivers efficient sterilization with a cycle time of approximately 30 minutes, solidifying its position as a leading choice in veterinary sterilization equipment.

[Read More](#)

SKU: AM-M11-022

Price: Request Quote



MIDMARK M9D MANUAL DOOR DENTAL AUTOCLAVE

The **Midmark M9D Manual Door Dental Autoclave (Refurbished)** is a reliable and efficient solution for dental instrument sterilization. Designed to ensure precision and safety, this autoclave features an intuitive start button that automatically fills the chamber with the correct amount of water, streamlining the sterilization process. It includes a manual door that can be opened to release steam, allowing for the fastest drying times and enhancing workflow efficiency. Built with safety in mind, it boasts critical features such as overheat protection, a pressure relief valve, and a flashing "Error" light to alert users of any issues. The self-program controls offer versatility by enabling users to customize cycle parameters to meet specific sterilization needs, while the green LCD display ensures easy monitoring of the autoclave's status. With overall dimensions of 20.13"D x 15.3"W x 15.8"H and a chamber capacity of 3.5 gallons, it accommodates a range of instruments. Weighing 70 lbs and requiring a 115V power supply, the M9D is both robust and user-friendly, making it an ideal choice for dental practices focused on reliability and performance.

[Read More](#)

SKU: AM-M9D-022

Price: Request Quote



TUTTNAUER 2540EK KWIKLAVE ELECTRONIC AUTOCLAVE (REFURBISHED)


The **Tuttnauer 2540EK Kwiklave Electronic Autoclave**, available in a refurbished model, is a versatile and efficient sterilization device ideal for various medical and laboratory settings. It offers four distinct programs—Unwrapped, Wrapped, Liquids, and a Drying cycle—facilitating comprehensive sterilization processes. Its intuitive one-touch operation allows electronics to manage the entire cycle seamlessly from start to finish, ensuring ease of use and reliability. The autoclave features a spacious 22-gallon chamber volume and a closed-door drying cycle, maintaining the sterility of instruments post-sterilization. Designed with user convenience in mind, it includes an automatic shut-off, audible and visual alerts for cycle interruptions, and a user-friendly keypad. For added utility, an optional printer is available for documenting sterilization cycles. All Tuttnauer autoclaves are A.S.M.E. stamped and U.L. approved, guaranteeing high standards of safety and performance. The autoclave has a chamber size of 10 inches in diameter and 19 inches in depth, accommodating four trays with each tray measuring 16.3 inches by 6.7 inches. Operating on 220-230V, this steam autoclave (Class S) is compact yet robust, with overall dimensions of 21.5 inches in depth, 20 inches in

width, and 14.4 inches in height, and has a right-side hinged door swing.

[Read More](#)

SKU: AM-2540EK

Price: Request Quote



CASSETTE FOR SCICAN STATIM G4 5000 AUTOCLAVE (REFURBISHED)


The **Cassette for Scican Statim G4 5000** is an essential accessory designed to enhance the efficiency of the Statim sterilization process. Engineered for optimal performance, this cassette accommodates a wide range of dental and surgical instruments, ensuring effective steam sterilization in a compact format. Its innovative design allows for rapid heat transfer, significantly reducing cycle times and improving workflow in busy practices. The cassette is constructed from durable materials that withstand repeated use, ensuring longevity and reliability. With easy loading and unloading features, it simplifies the sterilization process, making it an invaluable tool for healthcare professionals seeking to maintain high standards of hygiene and efficiency in their sterilization practices.

Grey Handle

[Read More](#)

SKU: 01-114989S-R

Price: Request Quote



CASSETTE FOR SCICAN STATIM G4 5000 AUTOCLAVE (NEW)

The **Cassette for Scican Statim G4 5000** is an essential accessory designed to enhance the efficiency of the Statim sterilization process. Engineered for optimal performance, this cassette accommodates a wide range of dental and surgical instruments, ensuring effective steam sterilization in a compact format. Its innovative design allows for rapid heat transfer, significantly reducing cycle times and improving workflow in busy practices. The cassette is constructed from durable materials that withstand repeated use, ensuring longevity and reliability. With easy loading and unloading features, it simplifies the sterilization process, making it an invaluable tool for healthcare professionals seeking to maintain high standards of hygiene and efficiency in their sterilization practices.

Grey Handle

[Read More](#)

SKU: 01-114989S

Price: Request Quote



STERI-DENT 200 DRY HEAT STERILIZER

The **Steri-Dent 200 Dry Heat Sterilizer** stands out as an exceptional method for sterilizing sharp instruments, effectively killing bacteria without moisture and preserving the sharpness of the tools. This versatile equipment not only serves as an efficient sterilizer but also functions as a drying unit for wet packs or as a reliable backup to an autoclave in a busy practice. One of its most appealing features is that it is maintenance-free, eliminating the need for routine cleaning that wet heat methods require. Safety is significantly enhanced as dry heat sterilization operates without the need for pressure, steam, or chemicals. The Steri-Dent 200 is equipped with a rapid-start heating element for quick operation, making it both safe and simple to use. Convenience is further underscored by utility trays that can be sterilized with lids on, ensuring they remain uncontaminated in storage. Automatic timers with a heat "hold" capability make it an ideal choice for continuous use. Constructed from high-quality stainless steel, the Steri-Dent 200 is built to last, resisting rust and ensuring long-term durability. Beyond its technical advantages, it is also value-priced, offering low purchase and operating costs, and eliminating the need for distilled water. The Steri-Dent 200 Dry Heat Sterilizer combines efficiency, safety, and economy, making it an indispensable tool for any medical practice.

[Read More](#)

SKU: AM-SD200

Price: Request Quote



TUTTNAUER T-TOP 10S GRAVITY AUTOCLAVE

The **Tuttnauer T-TOP 10S Gravity Autoclave** represents a significant advancement in sterilization technology, specifically designed for small clinics, laboratories, and other healthcare facilities. This model, known for its robust reliability and user-friendly interface, utilizes gravity displacement technology to effectively remove air and sterilize equipment using steam under pressure. Its compact design makes it ideal for environments where space is at a premium, yet it does not compromise on the capacity or the quality of sterilization. The T-TOP 10S also features a range of safety and efficiency protocols, including a heat-insulated double-locking door mechanism and an intuitive control panel that simplifies the sterilization process. This autoclave is an excellent choice for medical professionals seeking a cost-effective, efficient solution for ensuring the sterility of their instruments and supplies.

[Read More](#)

SKU: AM-T-TOP10

Price: Request Quote



TUTTNAUER LABSCI 85LV VERTICAL AUTOCLAVE

The **Tuttnauer LabSci 85LV Vertical Lab Autoclave** is a state-of-the-art steam sterilizer designed to meet the demanding needs of educational institutions and research facilities. With a generous 85-liter chamber volume, this top-loading autoclave maximizes the efficiency of sterilization processes by accommodating large and varied loads, including up to four 2000 mL Schott-Duran flasks or three 2000 mL Erlenmeyer flasks. It is equipped with four validated sterilization programs and 16 fully customizable options, making it suitable for sterilizing liquids, plastics, textiles, glass, and surgical instruments. The LabSci 85LV features a user-friendly, multi-colored touchscreen interface, making operation straightforward and intuitive. Additionally, its mobility is enhanced with swivel wheels, allowing easy relocation within the lab space without the need for plumbing. Its advanced features, such as integrated temperature probes, a built-in cooling fan, and memory for 999 cycles, along with compliance-enhancing connectivity options like USB, WiFi, and optional printer capabilities, ensure that this autoclave is not only practical but also



compliant with rigorous standards. The Tuttnauer LabSci 85LV is an essential tool for laboratories seeking a reliable, efficient, and economic sterilization solution.

[Read More](#)

SKU: AM-LABSCI85LV

Price: Request Quote



TUTTNAUER LABSCI 60LV VERTICAL AUTOCLAVE

The **Tuttnauer LabSci 60LV Vertical Lab Autoclave** represents a significant advancement in steam sterilization technology, particularly suited for educational institutions and research facilities. This top-loading autoclave, part of Tuttnauer's T-Lab ECO line, is engineered for economic and efficient sterilization of both liquids and solids, boasting a generous 60-liter chamber capacity. Its design facilitates ease of use with a height of just over three feet and a top-loading feature that eliminates the need for counter space. Enhanced mobility is provided by swivel wheels, allowing for easy relocation within the lab environment. The autoclave offers high functionality with four validated programs and sixteen customizable options, making it versatile for a variety of sterilization needs including plastics, textiles, glass, and instruments. Its user-friendly interface features a high-resolution touchscreen, and the device requires no complex plumbing for installation—simply plug it into a power source. With additional features like built-in temperature probes, a cooling fan, and memory for up to 999 cycles, the Tuttnauer LabSci 60LV is designed to optimize lab workflow and productivity, while ensuring rigorous sterilization standards are met.

[Read More](#)

SKU: AM-LABSCI60LV

Price: Request Quote





TUTTNAUER LABSCI15L+ BENCHTOP AUTOCLAVE

The **Tuttnauer LABSCI15L+ Benchtop Autoclave** stands as a premier solution for modern laboratories requiring reliable and efficient sterilization. With a generous 85L (approximately 17 gallons) chamber volume, this autoclave is equipped to handle a variety of materials—from glassware and instruments to biological wastes—across four standard program cycles and up to 16 customizable cycles. Crafted from robust 316L or 316Ti stainless steel, the chamber guarantees longevity and resistance to corrosion. The autoclave features an advanced multi-color control panel for ease of operation, and it requires no plumbing, simply plug in and start. Additional features such as a built-in printer provide detailed cycle documentation, enhancing traceability and compliance. The two-year warranty on parts and labor ensures dependable service while its array of features like fast cooling, efficient drying, and optional biohazard filters make it a versatile asset in any lab setting. Whether for research or clinical applications, the LABSCI 15L+ enhances productivity and maintains strict sterility standards with every use.

[Read More](#)

SKU: AM-LABSCI15L+

Price: Request Quote



TUTTNAUER LABSCI15L BENCHTOP AUTOCLAVE

The **Tuttnauer LABSCI 15L autoclave** represents a pinnacle of convenience and reliability for laboratory sterilization needs. Engineered to handle a broad array of materials - from liquids and media to instruments, glassware, clothing, and even waste - this versatile unit ensures a high standard of sterilization. Its "plug & play" design is particularly advantageous, requiring no complex plumbing installations, thus facilitating easy integration into any laboratory setting. With a focus on safety and efficiency, the LABSCI 15L provides rapid sterilization cycles, making it an indispensable tool for maintaining stringent hygiene and safety standards in scientific research and healthcare environments. Whether for routine or intensive sterilization tasks, the Tuttnauer LABSCI 15L is designed to deliver exceptional performance without compromise.

[Read More](#)

SKU: AM-LABSCI15L

Price: Request Quote



TUTTNAUER LABSCI11L BENCHTOP AUTOCLAVE

The **Tuttnauer LABSCI 11L autoclave** is an essential tool for laboratory environments, designed with versatility and convenience in mind. This innovative system effectively sterilizes everything required in a lab setting, including liquids, media, instruments, glassware, clothing, and even waste. Its plug-and-play design ensures easy installation and operation as it does not require any complex plumbing connections. With its reliable and rapid sterilization capabilities, the LABSCI 11L ensures that all materials are handled safely and are free from contaminants, thereby facilitating a sterile, efficient, and productive laboratory workspace.

[Read More](#)

SKU: AM-LABSCI11L

Price: Request Quote



TUTTNAUER DS1000 STEAM DISTILLER 1 GALLON

The **Tuttnauer DS1000 Steam Distiller 1Gal** stands as a premier choice for those in need of high-quality, distilled water in laboratory, medical, or industrial settings. This compact yet powerful distillation unit is designed to convert tap water into pure, contaminant-free steam distilled water at a capacity of 1 gallon per cycle. The DS1000 model is renowned for its reliability, efficiency, and ease of use, featuring an intuitive operation system that ensures consistent production of distilled water without the need for constant monitoring. Its robust construction and innovative design make it a durable option for various applications, providing users with a sustainable solution for their distilled water needs. Whether for preparing medical instruments, conducting laboratory experiments, or any application requiring the utmost purity in water, the Tuttnauer DS1000 Steam Distiller 1Gal delivers exceptional performance and quality assurance.

[Read More](#)

SKU: DS1000

Price: Request Quote



TUTTNAUER DS1G STEAM DISTILLER 4L - 1L PER HR

The **Tuttnauer DS1G Steam Distiller 4L** represents a pinnacle of purification technology, specifically designed for those who demand the utmost in water purity. This compact, yet powerful device is capable of producing 4 liters of distilled water in just a few hours, making it an ideal solution for laboratories, medical facilities, or even for personal use in homes where water quality is a paramount concern. Its steam distillation process is meticulously engineered to remove contaminants, bacteria, and other impurities, ensuring that the water it produces is not only safe but of the highest possible quality. With its robust build and efficiency, the Tuttnauer DS1G offers a reliable and convenient way to access distilled water, while its design and ease of use make it an attractive addition to any setting that values purity and health. Distills 1 Liter per hour.

[Read More](#)

SKU: DS1G

Price: Request Quote



TUTTNAUER 3870HSG-WS PRE/POST VACUUM 85L / 22.5 GALLON AUTOCLAVE

The **Tuttnauer 3870HSG-WS autoclave** is an 85 liter autoclave that is ideal for ASC's and large clinics. The 3870HSG-WS autoclave provides hospital grade sterilization and medium loading capacity at a cost-effective price that meets the highest international standards.

The 3870 HSG-WS is a hospital grade Class B autoclave, suitable for large medical & dental clinics, and operating rooms. Designed to be a fast cycle sterilizer that is mobile, operates 24 hours a day, and can optionally operate independent of building utilities.

Independent Mode - no need to connect the autoclave to the building utilities for distilled water, tap water or drain

[Read More](#)

SKU: 3870HSG-WS

Price: Request Quote



STERIDATA LOGGER STERILIZER PRINTER, AUTOCLAVE PRINTER, PAPERLESS - RITTER / MIDMARK M9 / M11 NEW SERIES

SteriData™ Logger for Midmark Ritter M9 and M11 Ultraclaves, M3 Ultrafast sterilizer, and SciCan Statim 2000 and 5000 cassette sterilizers*

The paperless alternative to a paper Printer

This SteriData Logger unit plugs in just like the paper printer and sits on top of your Sterilizer. Now all your sterilizer cycle data is saved to the SD card (included) automatically.

[Read More](#)

SKU: AM-M1219-DL

Price: \$1,095.00



TUTTNAUER 3870E AUTOCLAVE

The refurbished **Tuttnauer 3870E Autoclave** offers a spacious **84-liter (22-gallon) capacity within its 15" x 27" chamber**, accommodating substantial sterilization loads. Automated digital controls simplify operation, featuring a 60-minute timer and a clear digital display of temperature and pressure. The included support stand and two 14" x 26.1" stainless steel tray shelves enhance loading convenience. Safety features include a heat-insulated, double-locking door that prevents opening while the chamber is pressurized, protecting users from burns. An automatic low-water-level shut-off safeguards the heating elements from damage. With a temperature range of 105°C to 137°C (221°F to 279°F) and a heating time to 134°C (273°F) of 30 minutes, the 3870E ensures efficient sterilization. The autoclave is UL Listed, CSA Certified, and ASME Stamped, meeting stringent safety and quality standards, and comes with a 6-month parts warranty. Its overall compact dimensions of 34.5" x 26" x 20.6" make it suitable for various settings.

[Read More](#)

SKU: AM-3870E

Price: Request Quote



MIDMARK RITTER M11 AUTOCLAVE, M11-001

Refurbished Ritter M11 with 11" x 18" chamber. It's chamber size gives it the distinction of having the largest sterilizing capacity of any standard countertop sterilizer on the market, making it perfect for larger packs and instruments.

[Read More](#)

SKU: AM-M11-001

Price: Request Quote



TUTTNAUER 3870EA AUTOCLAVE GEN 2 (NEW)

Includes:
3870EA Autoclave
Two Instrument Trays and one holder
Autoclave Cleaner
Drain Tube
Operators Manual

Includes Printer

[Read More](#)

SKU: AM-3870EA-N

Price: Request Quote



3870EA AUTOCLAVE GEN 2 (REFURBISHED)

Includes:
3870EA Autoclave
Two Instrument Trays and one holder
Autoclave Cleaner
Drain Tube
Operators Manual

For Printer Request AM-TUT-P

[Read More](#)

SKU: AM-3870EA-G2

Price: Request Quote



**No Image
Available**

T-EDGE 10 OR 11 BARCODE PRINTER

The printer(s) are optional and can easily be installed and connected to the autoclave Prints a detailed history of each cycle performed. (This can be used for the record or for subsequent consideration.)

The printing is on thermal paper with a defined set of characters per line and contains important information such as some of the main following details:

- Date: Time: Ser. Num: Model: Version:
- Cycle Num: Cycle Name: Ster Temp: Ster Time: Dry Time: End Temperature

[Read More](#)

SKU: THE002-0116

Price: Request Quote



TUTTNAUER T-EDGE 11S GRAVITY AUTOCLAVE

Fast Cycles: 36 Minute wrapped instruments complete cycle, including drying

Optimal chamber capacity: More than 80% chamber capacity extra-long trays

Convenient loading: Open to the door to 110° ensures it will not block your smooth loading of the chamber

Remote Wi-Fi: Monitor you autoclave via smartphone or tablet app

Advanced documentation: A variety of documentation and traceability options and barcodes

User-friendly touch screen control panel: Displays programs, cycles times, water levels and more

[Read More](#)

SKU: AM-T-EDGE11S

Price: Request Quote



TUTTNAUER T-EDGE 11B PRE POST VACUUM AUTOCLAVE

Fast Cycles: 36 Minute wrapped instruments complete cycle, including drying

Optimal chamber capacity: More than 80% chamber capacity extra-long trays

Convenient loading: Open to the door to 110° ensures it will not block your smooth loading of the chamber

Remote Wi-Fi: Monitor you autoclave via smartphone or tablet app

Advanced documentation: A variety of documentation and traceability options and barcodes

User-friendly touch screen control panel: Displays programs, cycles times, water levels and more

[Read More](#)

SKU: AM-T-EDGE11B

Price: Request Quote





MIDMARK RITTER M11D DENTAL MANUAL AUTOCLAVE

The **Midmark Ritter M11D Dental Manual Autoclave** is a highly efficient and reliable sterilization solution, particularly noted for its refurbished model equipped with an expansive 11" x 18" chamber. This generous chamber size distinguishes it as the countertop sterilizer with the largest capacity currently available, making it exceptionally suitable for processing larger packs and instruments commonly used in dental practices. Key features of the Ritter M11 include a manual door opening designed specifically for dental applications and a robust Steam-Flush Pressure Pulse Air Removal System, which guarantees thorough steam penetration for optimal sterilization. The unit is equipped with self-program controls that allow users to tailor cycle parameters to specific instrument processing requirements, alongside pre-programmed controls for added convenience. The user-friendly green LCD display ensures easy visibility of all indicators, enhancing operational efficiency. The autoclave also includes two trays measuring 9" x 15" and two trays measuring 6.6" x 15", fitting within its external dimensions of 17.8" wide, 23.8" deep, and 17.8" high. Despite its substantial capacity, the Ritter M11 maintains a cycle time of about 30 minutes and a weight of 90 pounds, making it a practical addition to any dental practice seeking reliable and efficient sterilization solutions.

[Read More](#)

SKU: AM-M11D-022

Price: Request Quote



TUTTNAUER EZ 11 PLUS NO PRINTER

Closed door drying, meaning the door remains closed through the dry cycle to maintain sterility and ensure efficient drying. A high efficiency air pump circulates filtered air through the chamber for fast drying. 0.2µ HEPA filter (replaceable).

[Read More](#)

SKU: AM-EZ11PLUS

Price: Request Quote



TUTTNAUER EZ10K 220V AUTOCLAVE

The **Tuttnauer EZ10k 220V Autoclave** is a cutting-edge sterilization tool renowned for its efficiency and reliability in medical and laboratory settings. Designed with rapid cycle times, it features a fully automatic, user-friendly interface, making it ideal for high-demand environments where time is crucial. Its advanced closed-door drying system ensures instruments are thoroughly dried post-sterilization, preserving their integrity for immediate use. Constructed with durability in mind, the EZ10k includes robust safety features like automatic shut-off and pressure release mechanisms, which provide operator safety and consistent performance. The spacious chamber accommodates a variety of instrument sizes, allowing for diverse sterilization needs while upholding high infection control standards. The Kwiklave model, a quick cycle variant, reduces total sterilization time by 50%, making it perfect for quick instrument turnover or high-volume settings. Operating on a voltage between 220V and 235V, it may require a buck/boost transformer for optimal performance. With a chamber volume of 6 gallons and a power rating of 2200W, the EZ10k is both powerful and versatile, though it does not include a printer.

[Read More](#)

SKU: AM-EZ10K-R

Price: Request Quote



RITTER M11 VETERINARY AUTOCLAVE

Refurbished Ritter M11 with 11" x 18" chamber.

Its chamber size gives it the distinction of having the largest sterilizing capacity of any standard countertop sterilizer on the market, making it perfect for larger packs and instruments.

[Read More](#)

SKU: AM-M11-033

Price: Request Quote



TUTTNAUER ELARA 11 AUTOCLAVE

With over 80 years of expertise in sterilization and infection control, Tuttnauer has always provided reliable and advanced technology. The Elara 11, Tuttnauer's latest addition to the pre & post vacuum autoclave family, continues this tradition. It is an ideal sterilization solution for any dental, outpatient or private clinic which needs the highest sterilization standards.

The Elara 11 is designed for sterilization of wrapped, porous and hollow loads. The chamber volume is large, external dimensions are small, making it especially suitable for smaller clinics requiring a high turnover. The control system is smart and intuitive, easy to use and has a large memory for consistent and reliable performance.

[Read More](#)

SKU: AM-ELARA11-R

Price: Request Quote



RITTER M11 AUTOCLAVE, M11-001 - RENTAL

Refurbished Ritter M11 with 11" x 18" chamber. It's chamber size gives it the distinction of having the largest sterilizing capacity of any standard countertop sterilizer on the market, making it perfect for larger packs and instruments.

[Download Brochure](#)

[Read More](#)

SKU: N/A

Price: Request Quote

MIDMARK M9/M11 PRINTER AND CABLE

Midmark M9/M11 Printer and cable
Cable and 1 roll of recorder paper included

Compatible with models:
M9-020, M9-022, M9-033
M11-020, M11-022, M11-033

[Read More](#)

SKU: 9A259001

Price: Request Quote



RITTER MIDMARK M9, M11, OLD STYLE, CONNECTION CABLE INCLUDED (RETIRED)



*** This product is retired and no longer available: for a newer comparable product please see the [SteriData Logger Sterilizer Printer, Autoclave Printer, Paperless - Ritter / Midmark M9 / M11 NEW Series](#)

[Read More](#)

SKU: 9A155001

Price: Request Quote

STATIM 2000 G4 CASSETTE STERILIZER



The Statim 2000 G4 is the same renowned and trusted autoclave it has been for over 20 years but now boasts a new contemporary look and connectivity that is the first of its kind. The G4 technology will change the way you interact by providing a direct channel of communication through the Internet to you, or anyone you desire. Still powered by SciCan's signature steam technology to provide sterilization and dryness at speeds faster than conventional chambered autoclaves, the Statim has been drastically upgraded with a level of interactivity never seen before.

[Read More](#)

SKU: AM-STATIM-2000-G4

Price: Request Quote



TUTTNAUER 2540M AUTOCLAVE

The **Tuttnauer 2540M Autoclave** is a reliable, cost-effective sterilization solution designed for daily use with minimal maintenance. Featuring an automatic shut-off at the end of both sterilization and dry cycles, a dual safety thermostat for overheating protection, and a durable electro-polished chamber and door, this autoclave ensures safe and efficient operation. Its user-friendly design includes a front-mounted drain valve for quick reservoir drainage and a double safety locking device that prevents the door from opening while the chamber is pressurized. As one of our best-selling models, we keep the Tuttnauer 2540M in stock for same-day shipping, making it an ideal workhorse sterilizer for busy healthcare and laboratory environments.

[Read More](#)

SKU: AM-2540M

Price: Request Quote

TUTTNAUER 2540MK



The Tuttnauer M Series is also available in a quick cycle model, the Kwiklave (MK). This sterilizer offers all the proven, dependable features of our manual autoclave, with the added benefit of reducing the total heat up time. This autoclave is the perfect choice for quick instrument turn-around or a high volume office. Speed and reliability are accomplished without sacrificing load size. 1 year parts and labor warranty.

[Read More](#)

SKU: AM-2540MK

Price: Request Quote

COX 6000 RAPID DRY HEAT



The COX Rapid Heat Transfer Sterilizer is a fast, simple solution to your sterilization needs. Because it utilizes heated air that moves at a rate of 2,500 feet/minute, it is able to sterilize objects in as little as six (6) minutes -- faster and easier than any other sterilizer available.

No longer will you and your staff be bothered with all the problems inherent with other forms of sterilizations. Since the dry heat system won't corrode or rust your instruments, you won't have to replace them so often -- thus saving you money.

[Download Brochure](#)

[Read More](#)

SKU: AM-C6000

Price: Request Quote

DRI-CLAVE 3000 DRY HEAT STERILIZER 3 TRAY



The FDA approved technique of dry heat sterilization works without pressure, steam, or chemicals. Dry heat sterilizes in the shortest possible time and you never need to dry the instruments. Ideal for all medical and dental instruments, it sterilizes instruments, needles, syringes, glassware, etc. at 320F (160C). Sharp instruments remain sharp; no steam or moisture to corrode the finest edge. Both operate on 115V AC and are UL Listed and FDA Registered.

[Download Brochure](#)

[Read More](#)

SKU: AM-DC-EC-3000

Price: Request Quote



STATIM 5000 STERILIZER

The Statim 5000 sterilizer has an additional sterilization cycle and a roomier cassette designed to process bulkier instruments and greater loads. It's a reliable cassette sterilizer that's compact enough to fit in patient treatment areas and fast enough to sterilize expensive Lasik, PRK, endoscopic and cataract instruments, including phaco handpieces, between procedures - it can complete an unwrapped sterilization cycle in just nine minutes and a wrapped cycle in just 17.5 minutes.*

[Download Brochure](#)

[Read More](#)

SKU: AM-STATIM-5000-R

Price: Request Quote



TUTTNAUER 3870M AUTOCLAVE (REFURBISHED)

The **Tuttnauer 3870M Manual Autoclave** stands as the largest manual table-top autoclave offered by Tuttnauer, continuing to fulfill a significant demand in healthcare facilities and veterinary practices. This robust, straightforward unit provides essential sterilization capabilities, proving reliable even in challenging conditions. With a spacious 15" x 30" chamber, it efficiently accommodates longer instruments or a larger volume of smaller tools, making it an ideal choice for spay and neuter clinics, equine practices, bustling animal hospitals, and multi-practice medical offices. Designed as a freestanding unit, the 3870M requires a dedicated 230 Volt, 13.6 Amp power outlet and operates without the need for special plumbing, as its reservoir is filled with easily accessible distilled water. Routine maintenance is simplified with features like an automatic shut-off at the end of sterilization and dry cycles, a long-lasting electro-polished chamber and door, and a convenient front-located drain valve. Enhanced safety is ensured through a double safety locking device to prevent door opening under pressure and a dual safety thermostat safeguarding against overheating.

[Download Brochure](#)

[Read More](#)

SKU: AM-3870M

Price: Request Quote



SCICAN STATIM 5000 G4 CASSETTE AUTOCLAVE (REFURBISHED)

The SciCan Statim 5000 G4 Cassette Autoclave offers a fully automatic and exceptionally fast sterilization solution for busy practices. This refurbished unit utilizes four single-touch programs to optimize instrument cleaning from start to finish. With a larger cassette and an additional sterilization cycle, the Statim 5000 G4 can process larger instruments quickly and reliably, completing a wrapped cycle in just 17.5 minutes and an unwrapped cycle in a mere 9 minutes. Modernized with a user-friendly 3.5" color touchscreen and intuitive controls, the Statim 5000 G4 also features a built-in data logger, eliminating the need for messy thermal paper and external devices, saving cycle data to a USB drive or the cloud. With features like ultra-fast cycles, Dri-Tec Drying Technology, and aseptic transfer capabilities, the Statim 5000 G4 ensures exceptional instrument turnover, prolongs instrument life, and reduces instrument investment.

[Download Brochure](#)

[Read More](#)

SKU: AM-STATIM-5000-G4

Price: Request Quote



DRI-CLAVE DRY HEAT STERILIZER - 2 TRAY

The FDA approved technique of dry heat sterilization works without pressure, steam, or chemicals. Dry heat sterilizes in the shortest possible time and you never need to dry the instruments. Ideal for all medical and dental instruments, it sterilizes instruments, needles, syringes, glassware, etc. at 320 °F (160 °C). Sharp instruments remain sharp; no steam or moisture to corrode the finest edge. They operate on 115V AC and are UL Listed and FDA Registered. The best method for sharp, micro, delicate instruments: Dry heat sterilizer kills bacteria without moisture and will not rust or dull sharp or serrated instruments.


[Download Brochure](#)

[Read More](#)

SKU: AM-DC-EC-2000

Price: Request Quote





AMSCO EAGLE 10 PLUS

Amsco Eagle Ten Autoclave provides clinics with Safe, fast and dependable sterilizer. The Amsco Eagle Ten Autoclave is a large capacity, counter-top sterilizer perfect for any size practice. Amsco has been a world leader in sterilization equipment, and the Eagle Ten offers many sophisticated features. The spacious 10" x 16" chamber can hold two utensil trays per load for wrapped, unwrapped and liquid goods. The outside dimensions are compact, saving counter space in the work area. Simply select wrapped, unwrapped, packs or liquid sterilization with the touch of one button. Other features include optional programmable controls, self-diagnostic error codes and safe operation

[Download Brochure](#)

[Read More](#)

SKU: AM-AMSCO-EAGLE10+

Price: Request Quote



PELTON & CRANE VALIDATOR 10

The Pelton & Crane Validator 10 Autoclave is a economical, light-capacity sterilizing unit designed to accommodate the needs of small medical, dental, and other practices requiring reliable and safe access to sterilization equipment. Incorporated are a variety of safety and efficiency features, such as a door lock, which prevents the door from being opened while the chamber is pressurized. It also features a set of microprocessor-controlled features to ensure effective and safe operation; a vent valve will automatically open if chamber pressure exceeds specification with alarm and secondary release. Similar functionality exists for heat and electrical levels, as well as a sensor to prevent chamber overfilling. Unique "Pressure Locking" Door Design. Provides maximum safety.

[Download Brochure](#)

[Read More](#)

SKU: AM-PC-VALIDATOR-10

Price: Request Quote



MIDMARK RITTER M11 ULTRACLAVE (AUTOCLAVE), M11-022

The refurbished **Midmark Ritter M11 Ultraclave** (Autoclave) offers a compact yet spacious sterilization solution with its **11" x 18" chamber**, touted as the largest among standard countertop sterilizers. This increased capacity makes it ideal for processing larger instrument packs and equipment, improving workflow efficiency. The M11 features an automatic-opening door, facilitating ease of use and promoting rapid drying times, which is crucial for preventing contamination. Its Steam-Flush Pressure Pulse Air Removal System ensures thorough steam penetration for effective sterilization. The self-program controls allow customization of cycle parameters to meet specific instrument processing needs, offering flexibility for various materials and instruments. Pre-programmed controls provide convenient, standardized cycles for routine sterilization tasks. The green LCD display enhances visibility, allowing for easy monitoring of cycle indicators and parameters. The M11 Ultraclave combines capacity, convenience, and customizable control in a compact countertop design.

[Download Brochure](#)

[Read More](#)

SKU: AM-M11-022

Price: Request Quote

TUTTNAUER 2340M MANUAL AUTOCLAVE



The Tuttnauer manual autoclave is designed to compliment any healthcare facility. The bright new easy to use panel was designed with the operator in mind. Our 2340M autoclave will satisfy all of your sterilization needs without compromising on quality, safety or reliability.

[Download Brochure](#)

[Read More](#)

SKU: AM-2340M-R

Price: Request Quote



TUTTNAUER VALUEKLAVE 1730 AUTOCLAVE

The **Tuttnauer Valueklave 1730** will satisfy all of your basic sterilization needs at an economical price. The compact size (7" chamber) makes it a perfect primary unit for small offices or a second autoclave for larger practices. The fast cycle time and heated drying delivers a quick turnaround for wrapped or unwrapped instruments. The simple controls and easy closing door make this unit a breeze to operate. High speed and reliability are built into every Valueklave 1730.

[Download Brochure](#)

[Read More](#)

SKU: AM-1730

Price: Request Quote



PELTON & CRANE MAGNACLAVE

The Pelton & Crane Magnaclave is a fully self-contained steam sterilizer with a large chamber, making it ideal for busy group practice, emergency clinics, surgery centers and small hospitals. The chamber capacity allows it to hold up to 14 Hu-Friedy cassettes or equivalent. The unit is sold with one basket tray and one medium shallow tray. It is equipped with a safety valve which limits the maximum chamber pressure to 38 PSI. This unit is sold refurbished and includes a stand. Electrical Requirements: Dedicated 30 amp circuit/220-240 VAC or 208 VAC.

[Download Brochure](#)

[Read More](#)

SKU: AM-PC-MAGNACLAVE

Price: Request Quote

TUTTNAUER EZ10



The Tuttnauer EZ10 Fully Automatic Medical Autoclave Sterilizer features a simple one-touch design and pre-programmed cycles for: unwrapped instruments, wrapped and packs, glassware and an additional drying cycle. The Tuttnauer EZ10 Fully Automatic Autoclaves fill, sterilize, exhaust and dry at the touch of a button. Shut off at the end of both the sterilization and drying cycles is automatic. In addition, all program parameters can be customized and stored. The Tuttnauer EZ10 also features a replaceable HEPA filter air filter, and a high efficiency air pump that circulates hot air after the sterilization cycle for fast drying.

[Download Brochure](#)

[Read More](#)

SKU: AM-EZ10

Price: Request Quote



RITTER MIDMARK M9, M9-001

The M9 autoclave couldn't be easier and is time tested for reliability. You simply load it - set it - and forget it. Select your cycle and the time is set for you. Press "Start" and the M9 fills with the correct amount of water. Once the sterilization process is complete, the door opens automatically and quietly to dissipate the steam, giving you fast drying time for your instruments.

[Download Brochure](#)

[Read More](#)

SKU: AM-M9-001

Price: Request Quote

TUTTNAUER 3870EA AUTOCLAVE



The Tuttnauer 3870EA autoclaves have a closed-door active drying system to maintain sterility and ensure efficient drying of packs and pouches. A high efficiency air pump circulates hot air after the sterilization cycle for fast drying. These steam sterilizers use a replaceable HEPA 0.2m air filter.


[Download Brochure](#)

[Read More](#)

SKU: AM-3870EA

Price: Request Quote

STATIM 2000 STERILIZER



The STATIM 2000 Classic boasts a 6-minute cycle, but has the added advantage of employing a fully removable chamber so you can aseptically transport a set of freshly sterilized instruments directly to the procedure area.

The STATIM 2000 is fully automatic. Simply select one of the three sterilization programs and with a single touch of the keypad, begin the cycle. An advanced microprocessor controls all functions and maintains optimal sterilization conditions, from start to finish.

The STATIM 2000 is a reliable cassette sterilizer that's compact enough to fit in procedure areas and fast enough to sterilize jewelry, expensive Lasik, PRK, endoscopic and cataract instruments, including phaco hand pieces, between procedures - it can complete an unwrapped sterilization cycle in just six minutes and a wrapped cycle in 14 minutes.

[Download Brochure](#)

[Read More](#)

SKU: AM-STATIM-2000-R

Price: Request Quote

RITTER M11 AUTOCLAVE, M11-001

Refurbished Ritter M11 with 11" x 18" chamber. It's chamber size gives it the distinction of having the largest sterilizing capacity of any standard countertop sterilizer on the market, making it perfect for larger packs and instruments.

[Download Brochure](#)

[Read More](#)

SKU: AM-M11

Price: Request Quote



TUTTNAUER 2340EA AUTOCLAVE

This table-top autoclave is designed for the sterilization of wrapped and unwrapped instruments, glassware and other related items found in dental, medical and veterinary clinics, tattoo parlors, first aid rooms, hospitals, laboratories and other related industries.

[Download Brochure](#)

[Read More](#)

SKU: AM-2340EA-R

Price: Request Quote



TUTTNAUER EZ9 PLUS



The newly designed, ergonomically friendly EZGlide door features a sleek modern look and effortless opening and closing.

Innovative EZView Multi-Color Display offers real time cycle information that is easy to read, even from a distance. The EZView displays cycle parameters, temperature, and pressure in the chamber. During a cycle each stage is displayed with an active color wheel and text.

Closed door drying, meaning the door remains closed through the dry cycle to maintain sterility and ensure efficient drying. A high efficiency air pump circulates filtered air through the chamber for fast drying. 0.2 μ HEPA filter (replaceable).

[Read More](#)

SKU: N/A

Price: Request Quote

INDEX

3

3870EA Autoclave Gen 2 (Refurbished) 20

A

Alternative Data Logger to Midmark 9A682001 for M9/M11-001 through - 019 5

Alternative Data Logger to Midmark 9A682001 for M9/M11-020 through - 022 5

Alternative Data Logger to Midmark 9A682001 for M9/M11-040 through - 042 4

Amsco Eagle 10 Plus 31

C

Cassette for Scican Statim G4 5000 Autoclave (New) 12

Cassette for Scican Statim G4 5000 Autoclave (Refurbished) 11

Cox 6000 Rapid Dry Heat 28

D

Dri-Clave 3000 Dry Heat Sterilizer 3 Tray 28

Dri-Clave Dry Heat Sterilizer - 2 Tray 30

E

Enbio PRO Class B Autoclave 6

Enbio S Class B Autoclave 6

M

Midmark M9/M11 Printer and cable 25

Midmark M9D Manual Door Dental Autoclave 10

Midmark Ritter M11 Autoclave, M11-001 19

Midmark Ritter M11 Ultraclave (Autoclave), M11-022 32

Midmark Ritter M11 UltraClave Veterinary Autoclave, M11-033 9

Midmark Ritter M11D Dental Manual Autoclave 22

Midmark Ritter M9 UltraClave (Autoclave), M9-022 8

P

Pelton & Crane Magnaclave 33

Pelton & Crane Validator 10 31

R

Ritter M11 Autoclave, M11-001 36

Ritter M11 Autoclave, M11-001 - RENTAL 25

Ritter M11 Veterinary Autoclave 24

Ritter Midmark M9, M11, Old Style, Connection Cable Included (Retired) 26

Ritter Midmark M9, M9-001 34

S

Scican Statim 5000 G4 Cassette Autoclave (New) 8

Scican Statim 5000 G4 Cassette Autoclave (Refurbished) 30

Statim 2000 G4 Cassette Sterilizer 26

Statim 2000 Sterilizer 35

Statim 5000 Sterilizer 29

Steri-Dent 200 Dry Heat Sterilizer 12

SteriData Logger Paperless Sterilizer Cycle Printer For Statim 2
SteriData Logger Sterilizer Printer, Autoclave Printer, Paperless - Ritter / Midmark M9 / M11 Gen1
Series 2
SteriData Logger Sterilizer Printer, Autoclave Printer, Paperless - Ritter / Midmark M9 / M11 NEW Series
18

T

T-Edge 10 or 11 Barcode Printer 20
Tuttnauer 2340EA Autoclave 36
Tuttnauer 2340M Manual Autoclave 32
Tuttnauer 2540M Autoclave 27
Tuttnauer 2540MK 27
Tuttnauer 3870E Autoclave 18
Tuttnauer 3870EA Autoclave 35
Tuttnauer 3870EA Autoclave Gen 2 (New) 19
Tuttnauer 3870HSG-WS Pre/Post Vacuum 85L / 22.5 Gallon Autoclave 17
Tuttnauer 3870M Autoclave (New) 7
Tuttnauer 3870M Autoclave (Refurbished) 29
Tuttnauer DS1000 Steam Distiller 1 Gallon 16
Tuttnauer DS1G Steam Distiller 4L - 1L per Hr 17
Tuttnauer Elara 11 Autoclave 24
Tuttnauer EZ 11 PLUS No Printer 23
Tuttnauer EZ10 34
Tuttnauer EZ10k 220V Autoclave 23
Tuttnauer EZ9 Plus 37
Tuttnauer LABSCI 60LV Vertical Autoclave 14
Tuttnauer LABSCI 85LV Vertical Autoclave 13
Tuttnauer LABSCI11L Benchtop Autoclave 16
Tuttnauer LABSCI15L Benchtop Autoclave 15
Tuttnauer LABSCI15L+ Benchtop Autoclave 15
Tuttnauer T-Edge 11B Pre Post Vacuum Autoclave 21
Tuttnauer T-Edge 11S Gravity Autoclave 21
Tuttnauer T-TOP 10B Class B Pre/Post Vac Autoclave 3
Tuttnauer T-TOP 10S Gravity Autoclave 13
Tuttnauer T-TOP 11B SFPP Autoclave 4
Tuttnauer T-TOP 11S SFPP Autoclave 7
Tuttnauer Valueklave 1730 Autoclave 33
Tuttnauer 2540EK Kwiklave Electronic Autoclave (Refurbished) 10