

TABLE OF CONTENTS

Bladder Scanners	2
Roll Stand for Infinium Vioscan Bladder Scanner	3
Infinium Vioscan Bladder Scanner	3
Verathon BVM 9500 BladderScan	4
Verathon BladderScan BVI 9600	4
Verathon BladderScan BVI 9400	5



BLADDER SCANNERS

Auxo Medical has bladder scanners that are portable, hand-held ultrasound device, which can perform a quick, easy and non-invasive scan of the bladder. With the bladder scanner you save time since the complete scanning lasts only a few minutes. It also does not require the intervention of an ultrasound technician and can prevent unnecessary invasive procedures. We provide trustworthy and safe brands like Verathon and Bladderscan.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[**BRAND NEW**](#)

[**RECONDITIONED**](#)

[**REFURBISHED**](#)



ROLL STAND FOR INFINIUM VIOSCAN BLADDER SCANNER

The **Infinium Vioscan Bladder Scanner Roll Stand** is designed to provide a mobile and convenient solution for healthcare professionals. This sturdy and reliable roll stand features a locking mechanism on the wheels to ensure stability during use, preventing any unwanted movement. The integrated probe holder ensures that the ultrasound probe is always within easy reach and securely stored when not in use. The additional storage shelf offers ample space for essential accessories and supplies, making it easier to keep everything organized and accessible during patient assessments.

The Infinium Vioscan Bladder Scanner Roll Stand is ideal for use in a variety of healthcare settings, including hospitals, clinics, and outpatient facilities. Its mobile design allows it to be easily transported between patient rooms or different departments, ensuring that the bladder scanner is always available where it is needed most.

[Read More](#)

SKU: 267.906.0001



INFINIUM VIOSCAN BLADDER SCANNER

The **Infinium VioScan™** is a state-of-the-art, portable, non-invasive ultrasound bladder scanning device designed to deliver precise bladder volume measurements. With its advanced algorithms, the VioScan ensures accurate assessments and is equipped with a built-in strip recorder for easy record-keeping. The device features an 8.4-inch high-definition color screen and a simplified, user-friendly interface, making it accessible for medical professionals. In expert mode, the VioScan offers real-time HD bladder images, while express mode allows for quick and efficient bladder urine volume measurement. Fast scanning is facilitated by built-in bladder border recognition, and auto alignment provides visual guides for real-time pre-scan feedback to ensure proper probe placement. The scanner also boasts a large, upgradable SD data storage capacity, capable of holding over 1000 cases, and includes settings for male, female, and pediatric patients. It supports multiple languages, including English, French, German, Spanish, Portuguese, Swedish, Norwegian, and Danish. Additional features include patient bar code scan capability, two rechargeable battery sets, and an optional rolling cart for enhanced mobility.

Comparable to Verathon BVI Bladder Scanner and BioCon® Bladder Scanner

[Read More](#)

SKU: 000.901.0001



VERATHON BVM 9500 BLADDERSCAN

The **Verathon BVM 9500 BladderScan** is a state-of-the-art ultrasound device designed for non-invasive bladder volume assessment. This portable and user-friendly system is essential for healthcare professionals in various settings, including hospitals, clinics, and long-term care facilities. The BVM 9500 enhances patient care by providing accurate and timely bladder measurements, helping to prevent complications associated with bladder overdistension and urinary retention.

[Read More](#)

SKU: AM-BVM9500



VERATHON BLADDERSCAN BVI 9600

The **Verathon BladderScan BVI 9600** is an innovative, portable ultrasound system designed to deliver quick and precise bladder volume measurements in a non-invasive manner. This device significantly reduces the risk of catheter-associated urinary tract infections (CAUTIs) by offering a safe alternative to traditional catheterization, thereby enhancing patient comfort. Its user-friendly design includes a single-button operation, making it an ideal choice for diverse healthcare settings. Key features of the BladderScan BVI 9600 include its ability to provide non-invasive bladder volume measurements, a real-time bladder volume display for immediate results, and a clear, easy-to-read display that ensures accuracy in interpreting results. The device's portability and lightweight construction allow for easy maneuverability, while its durable build ensures it can withstand the rigors of daily use. Minimal training is required due to its intuitive user interface, which boosts efficiency by delivering quick, accurate results. Overall, the Verathon BladderScan BVI 9600 is a cost-effective solution that not only reduces CAUTI risks but also improves patient comfort and streamlines healthcare workflows,

ultimately reducing costs associated with catheterization and its potential complications.

[Read More](#)

SKU: AM-BVM9600



VERATHON BLADDERSCAN BVI 9400

The **Verathon BladderScan BVI 9400** is a cutting-edge medical device renowned for its efficiency and ease of use in measuring bladder volume. With just a press of the scan button, this device swiftly captures ultrasonic reflections across multiple planes within the body, creating a three-dimensional image that accurately calculates and displays bladder volume—all without the need for a sonographer. Its versatility is evident in its distinct scan modes tailored for men, women, and small children, catering to diverse patient needs. The handheld probe is ergonomically designed for precision and comfort, while the compact console features a tilting color LCD screen and onboard printer, enhancing usability. The BVI 9400 integrates seamlessly with office systems via HIPAA-compliant ScanPoint® technology, allowing for easy data transfer, viewing, and archiving. Furthermore, it offers valuable tools such as a savings calculator, voice annotation, and a four-minute video tutorial that simplifies training. With its durable design, optional mobile cart, and wireless capabilities for EHR integration, the BladderScan BVI 9400 stands out as a reliable and indispensable tool in modern medical practice.

[Download Brochure](#)

[Read More](#)

SKU: AM-BVI-9400

INDEX

I

Infinium Vioscan Bladder Scanner 2

R

Roll Stand for Infinium Vioscan Bladder Scanner 2

V

Verathon BladderScan BVI 9400 4

Verathon BladderScan BVI 9600 3

Verathon BVM 9500 BladderScan 3