

TABLE OF CONTENTS

Carts / Transport	2
Harloff M-Series Emergency Cart - MDS3024B14-EMG1Q	2
Harloff M-Series Anesthesia Cart - MDS3030E15	3
4 Shelf Utility Cart 24"x48"x69"	3
3 Shelf Utility Cart 18"x36"x33"	4
Small Distribution Cart - OPEN CASE CART	4
Harloff M-Series Cart - Standard Width Short Cabinet with 24" of Vertical Drawer Space, Breakaway Lock, 4 Drawers 2-3, 1-6, 1-12 with Standard Package	5
Pedigo Roll-out wire shelf, with (2) position spring loaded shelf lock for Case Cart Models RCC-242 and RCC-245	5
Harloff A-Series Aluminum Cart - Standard Width Tall Cabinet with 30" of Vertical Drawer Space, Basic Electronic Pushbutton Lock with Key Lock Override	6
Large Distribution Cart - OPEN CASE CART	6
Medium Distribution Cart - OPEN CASE CART	7
5 Drawer Emergency Cart with bottom Panel	7
Case Cart, Closed, 28.125inD X 46inW X 39.375inH 2 Door Clean indicator	8
Mobile Base Cabinet	8
Quick Ship ICU 4 Drawer Iso Cart w/ Key Lock	9



CARTS / TRANSPORT

Finding the right utility cart for your facilities operational and clinical support need is important to having the right supplies at the bed side or transported to your sterile processing department. Auxo Medical offers a variety of Emergency Crash Carts, Anesthesia Carts, Malignant Hypothermia Carts, Narcotics Cabinets, Stainless Case Carts, Open Case Carts, and many more.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[**BRAND NEW**](#)

[**RECONDITIONED**](#)

[**REFURBISHED**](#)



HARLOFF M-SERIES EMERGENCY CART - MDS3024B14-EMG1Q

The **Harloff M-Series Emergency Cart (MDS3024B14-EMG1Q)** is a crucial asset for any medical facility needing a rapid response solution. This standard-width, short cabinet cart provides 24 inches of vertical drawer space, configured with two 3-inch drawers, one 6-inch drawer, and one 12-inch drawer, all secured by a breakaway lock. The included Emergency Package equips the cart with an IV pole, a side-mounted drop shelf, two utility hooks, a cardiac board, an oxygen tank holder, and two sets of adjustable drawer dividers. Standard features such as an antimicrobial plastic top, ergonomic push handles, smooth-rolling casters, full-extension drawer slides, protective bumpers, and breakaway lock seals ensure functionality and security. The cart is finished in a highly visible red paint.

[Read More](#)

SKU: MDS3024B14-EMG1Q

Price: Request Quote



HARLOFF M-SERIES ANESTHESIA CART - MDS3030E15

The **Harloff M-Series Anesthesia Cart, model MDS3030E15**, is designed for efficient and organized storage of anesthesia supplies. It features a standard configuration of five drawers: two 3-inch drawers, two 6-inch drawers, and one 12-inch drawer. The package includes an antimicrobial plastic top for a hygienic work surface and ergonomic push handles on one side for easy maneuverability. The cart is equipped with 5-inch ball-bearing Tente casters for smooth rolling and full-extension drawer slides for complete access to drawer contents. Drawer stabilizing pontoon bumpers ensure quiet and secure drawer operation. This model replaces the 6415E anesthesia cart.

[Read More](#)

SKU: MDS3030E15

Price: Request Quote



4 SHELF UTILITY CART 24"X48"X69"

Round Post 4-Shelf Wire Shelving Carts, 24"D x 48"W x 63"H; Cart kit includes (4) posts, (4) shelves, (16) collets, (2) 5" swivel casters and (2) 5" swivel casters with brakes; High quality, plated wire shelves and corner posts; Shelves are 24"D x 48"W and feature truss reinforced edges; Posts feature groove rings in 1" increments for shelf height adjustability; 1200 lb. total cart capacity

[Read More](#)

SKU: AM-R244863CC-4

Price: Request Quote



3 SHELF UTILITY CART 18"X36"X33"

Round Post 3-Shelf Wire Shelving Carts with shelf size 18"D x 36"W; Cart kit includes (3) shelves, (2) U-handles, (12) collets, (2) 5" swivel casters and (2) 5" swivel casters with brakes; High quality, plated wire shelves and corner posts with integrated handles; Shelves feature truss reinforced edges; Posts feature groove rings in 1" increments for shelf height adjustability; 1200 lb. total cart capacity

[Read More](#)

SKU: AM-R1836CUC-3

Price: Request Quote



Note: Image shown is representative of the general construction of this model. Please refer to specifications for accurate model attributes. Shown with optional bins.

SMALL DISTRIBUTION CART - OPEN CASE CART

38 1/2 W x 28 3/4 D x 38 1/2 H - Accommodates up to 5 shelves

[Read More](#)

SKU: AM-7910

Price: Request Quote



HARLOFF M-SERIES CART - STANDARD WIDTH SHORT CABINET WITH 24" OF VERTICAL DRAWER SPACE, BREAKAWAY LOCK, 4 DRAWERS 2-3, 1-6, 1-12 WITH STANDARD PACKAGE

M-Series Cart - Standard Width Short Cabinet with 24" of Vertical Drawer Space, Breakaway Lock, 4 Drawers 2-3, 1-6, 1-12. Standard Package Includes: Antimicrobial Plastic Top, Ergonomic Push Handles (one side), 5" Ball- Bearing Tente Casters, Full-Extension Drawer Slides, Drawer Pulls, Protective Stabilizing Pontoon Bumpers and 100 Breakaway Lock Seals. REPLACES 6300

[Read More](#)

SKU: MDS3024B14

Price: Request Quote



PEDIGO ROLL-OUT WIRE SHELF, WITH (2) POSITION SPRING LOADED SHELF LOCK FOR CASE CART MODELS RCC-242 AND RCC-245

Roll-out wire shelf, with (2) position spring loaded shelf lock. Shelf locks in storage and fully extended position, for RCC-242 and RCC-245.
Est. Lead Time (Business Days)

[Read More](#)

SKU: RCC-242-WSRO

Price: Request Quote



HARLOFF A-SERIES ALUMINUM CART - STANDARD WIDTH TALL CABINET WITH 30" OF VERTICAL DRAWER SPACE, BASIC ELECTRONIC PUSHBUTTON LOCK WITH KEY LOCK OVERRIDE

A-Series Aluminum Cart - Standard Width Tall Cabinet with 30" of Vertical Drawer Space, Basic Electronic Pushbutton Lock with Key Lock Override, Factory Set User and Supervisor Code, Cabinet May Be Programmed in the Field. 6 Drawers 3-3, 2-6, 1-9. Standard Package Includes: Antimicrobial Plastic Top, Ergonomic Push Handles (one side), 5" Ball-Bearing Tente Casters, Soft-Close Full Extension Drawer Slides, Drawer Pulls, and Protective Stabilizing Pontoon Bumpers. Replaces Harloff AL810E6

[Read More](#)

SKU: MPA3030E06

Price: Request Quote



Note: Image shown is representative of the general construction of this model. Please refer to specifications for accurate model attributes. Shown with optional bins.

LARGE DISTRIBUTION CART - OPEN CASE CART

32 1/2 W x 32 3/4 D x 61 1/2 H - Accommodates up to 9 shelves
6 in Casters, All swivel, 2 w/ Brakes

[Read More](#)

SKU: AM-7930

Price: Request Quote



Note: Image shown is representative of the general construction of this model. Please refer to specifications for accurate model attributes. Shown with optional bins.

MEDIUM DISTRIBUTION CART - OPEN CASE CART

42 1/2 W x 30 3/4 D x 44 H - Accommodates up to 6 shelves
6 in Casters, All swivel, 2 w/ Brakes

[Read More](#)

SKU: AM-7920

Price: Request Quote

5 DRAWER EMERGENCY CART WITH BOTTOM PANEL



Emergency Carts - ER Crash Carts
DRAWER EMERGENCY CART WITH Panel as bottom
Drawer config. 4-3"1-6"1-9"panel
27" DRAWER SPACE
42-3/4"H x 25"D x 32"W

[Read More](#)

SKU: AM-MBS-527P

Price: Request Quote

CASE CART, CLOSED, 28.125IND X 46INW X 39.375INH 2 DOOR CLEAN INDICATOR



Surgical Case cart, 28.125inD X 46inW X 39.375inH / 25inD
Internal: 40.25inW X 30.25inH

Download Brochure 

[Read More](#)

SKU: MMCC3-1PS-SSD-CDI

Price: Request Quote

MOBILE BASE CABINET



20D x 28W x 34H

Marine edge

(2) 4 Swivel & (2) 4 Locking Casters

[Read More](#)

SKU: WBC1824

Price: Request Quote

QUICK SHIP ICU 4 DRAWER ISO CART W/ KEY LOCK



New 4 drawer anesthesia cart with key lock

Drawer Configuration: 1-3in, 3-9in

[Read More](#)

SKU: AM-MKT-430-Y

Price: Request Quote

INDEX

3

3 Shelf Utility Cart 18"x36"x33" 3

4

4 Shelf Utility Cart 24"x48"x69" 3

5

5 Drawer Emergency Cart with bottom Panel 7

C

Case Cart, Closed, 28.125inD X 46inW X 39.375inH 2 Door Clean indicator 7

H

Harloff A-Series Aluminum Cart - Standard Width Tall Cabinet with 30" of Vertical Drawer Space, Basic Electronic Pushbutton Lock with Key Lock Override 5

Harloff M-Series Anesthesia Cart - MDS3030E15 2

Harloff M-Series Cart - Standard Width Short Cabinet with 24" of Vertical Drawer Space, Breakaway Lock, 4 Drawers 2-3, 1-6, 1-12 with Standard Package 4

Harloff M-Series Emergency Cart - MDS3024B14-EMG1Q 2

L

Large Distribution Cart - OPEN CASE CART 6

M

Medium Distribution Cart - OPEN CASE CART 6

Mobile Base Cabinet 8

P

Pedigo Roll-out wire shelf, with (2) position spring loaded shelf lock for Case Cart Models RCC-242 and RCC-245 5

Q

Quick Ship ICU 4 Drawer Iso Cart w/ Key Lock 8

S

Small Distribution Cart - OPEN CASE CART 4