

# TABLE OF CONTENTS

<b>Centrifuges</b> .....	3
Drucker Diagnostics DASH Coag 2-3-5 Minute STAT Coag Centrifuge .....	4
Drucker Diagnostics SERO 12 Blood Banking Centrifuge .....	4
Drucker Diagnostics HORIZON 24 Flex Programmable Routine Centrifuge .....	5
Drucker Diagnostics HORIZON 12 Flex Programmable Routine Centrifuge .....	6
Drucker Diagnostics TrueBond Veterinary centrifuge for blood, urine, fecal, and microhematocrit .....	6
Drucker Diagnostics HORIZON 6 Flex FA Programmable Routine Centrifuge .....	7
Drucker Diagnostics HORIZON 6 Flex Programmable Routine Centrifuge .....	8
Drucker Diagnostics HORIZON 24 High-Capacity Routine Centrifuge .....	8
Drucker Diagnostics HORIZON 12 Compact Routine Centrifuge .....	9
Drucker Diagnostics TrueBond Set 3 Set-and-lock veterinary centrifuge for blood, urine, fecal, and microhematocrit .....	10
Drucker Diagnostics HORIZON 6 Compact Routine Centrifuge .....	10
Drucker Diagnostics DASH Flex 24 Programmable STAT Centrifuge .....	11
Drucker Diagnostics DASH Flex 12 Programmable STAT Centrifuge .....	12
Drucker Diagnostics DASH Flex 6 DC Programmable STAT Centrifuge .....	12
Drucker Diagnostics DASH Apex 24 High Volume STAT Centrifuge .....	13
Drucker Diagnostics DASH Apex 12 Set-and-Lock STAT Centrifuge .....	14
Drucker Diagnostics 642M 12V DC Mobile Centrifuge .....	14
Drucker Diagnostics DASH Apex 6 Compact STAT Centrifuge .....	15
Drucker Diagnostics TrueBond Duet Compact Dual Rotor Veterinary Centrifuge .....	16
Drucker Diagnostics DASH Apex 4 Compact STAT Centrifuge .....	16
Drucker Diagnostics 642E Single-Speed Electronic Centrifuge .....	17
Drucker Diagnostics 642B Single-Speed Centrifuge .....	18
Drucker Diagnostics 614B Fixed Angle Single-Speed Centrifuge .....	18
StatSpin Express 4 Centrifuge .....	19
Serofuge 2002 Centrifuge .....	20





## CENTRIFUGES

**Auxo Medical offers a curated selection of centrifuges designed for various laboratory and medical applications. Our inventory features the StatSpin Express 4 and the Clay Adams Sero-Fuge 2002. These centrifuges are essential tools for separating components of liquids, such as blood, urine, or other biological samples, enabling further analysis and diagnosis.**

The StatSpin Express 4 is a high-speed horizontal centrifuge known for its rapid processing capabilities, enhancing laboratory efficiency. It is ideal for applications requiring quick separation of samples, such as blood cell separation for hematology analysis. The

Clay Adams Sero-Fuge 2002, on the other hand, is a versatile centrifuge designed for a broader range of applications, including serological testing, blood banking, and general laboratory procedures.

Centrifuges offered through Auxo Medical are known for their reliability, ease of use, and safety features. They are designed to meet the demands of busy laboratories, providing accurate and consistent results. The benefits of using these centrifuges include time-saving efficiency, enhanced sample quality, and improved diagnostic accuracy. Order your centrifuges through Auxo Medical today!

**Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:**

**[BRAND NEW](#)**

[RECONDITIONED](#)

[REFURBISHED](#)



## DRUCKER DIAGNOSTICS DASH COAG 2-3-5 MINUTE STAT COAG CENTRIFUGE

The **Drucker Diagnostics DASH Coag 2-3-5 Minute STAT Coag Centrifuge** is a cutting-edge device designed to optimize the efficiency and accuracy of coagulation testing. Engineered for rapid processing, this centrifuge is ideal for STAT labs where time is of the essence. It can process samples in as little as two minutes, ensuring that critical results are available promptly. The DASH Coag Centrifuge is equipped with intuitive digital controls, allowing technicians to easily adjust and monitor cycle parameters. Its robust design accommodates various tube sizes, making it versatile for different laboratory needs. With features like an automatic lid release and programmable cycles, it enhances user convenience and reduces the potential for handling errors. Built with precision in the USA, this centrifuge supports high-throughput environments while maintaining reliability and consistency in test outcomes. With its compact dimensions and lightweight design, it fits seamlessly into any laboratory setting, offering unparalleled performance without compromising on space.

[Read More](#)

**SKU:** 00-676-009-001

**Price:** \$3,695.00



## DRUCKER DIAGNOSTICS SERO 12 BLOOD BANKING CENTRIFUGE

The **Drucker Diagnostics SERO 12 Blood Banking Centrifuge** is an advanced and reliable solution for laboratories focused on high-quality blood analysis procedures. Engineered to accommodate a variety of tube sizes, the SERO 12 seamlessly performs blood grouping, typing, cross matching, and cell-washing tasks, ensuring efficiency and precision. Leveraging its heritage as an evolution of the world's most popular blood banking centrifuge, it processes samples swiftly without requiring individual handling. The user-friendly digital display allows for quick adjustments, offering up to 10 programmable cycles to tailor centrifuge operations to specific needs. Compatible with all Sero-Fuge™ rotors, the SERO 12 comes equipped with a 12-place swing-out rotor, providing exceptional value. Key features include a digital timer, customizable speed and time settings, an auto-pop lid for easy access, and a QR code for quick access to support materials. With a capacity for 6 capped or 12 uncapped tubes and a maximum speed of 3600 RPM, the SERO 12 is robust yet compact, weighing just 27 lbs. Manufactured in the USA and supported by a lifetime warranty on its rotor, the SERO 12 promises durability and dependable performance for demanding blood banking applications.

[Read More](#)

**SKU:** 00-483-009-001

**Price:** \$4,850.00



## DRUCKER DIAGNOSTICS HORIZON 24 FLEX PROGRAMMABLE ROUTINE CENTRIFUGE

The **HORIZON 24 Flex Programmable Routine Centrifuge by Drucker Diagnostics** is an indispensable asset for any laboratory, offering the versatility and reliability needed for demanding workloads. Engineered to handle all blood and urine tubes up to 50 mL effortlessly, this centrifuge stands out with its easy-to-use digital display, allowing users to program and store up to ten custom settings or select from three expertly designed presets. Its intuitive LED lid lighting system ensures that samples are never forgotten, enhancing sample quality by clearly indicating the centrifuge's operational status. The centrifuge's patented 2-stage horizontal rotor technology ensures optimal sample separation, while its convection-based cooling system maintains a consistently cool operation even during intensive use. Safety is paramount, with features such as an automatic lid lock, a clear shatter-proof lid for safe observation, and a reinforced guard bowl for enhanced durability. Powered by a brushless DC motor, this centrifuge promises quiet and maintenance-free operation, adhering to rigorous safety standards across

the US, Canada, and Europe. Compact yet robust, the HORIZON 24 Flex delivers efficient performance with a sound level of just 64 dBA, and is backed by a lifetime warranty on the rotor, making it a reliable choice for high-volume laboratory environments.

[Read More](#)

**SKU:** 00-384-009-000

**Price:** \$3,995.00



## DRUCKER DIAGNOSTICS HORIZON 12 FLEX PROGRAMMABLE ROUTINE CENTRIFUGE

The **Drucker Diagnostics HORIZON 12 Flex Programmable Routine Centrifuge** is a versatile and compact solution designed to meet the varied processing needs of routine laboratories, particularly where space is limited. This centrifuge allows users to program and store up to ten custom settings for time, RPM, and g-force, offering flexibility and convenience through its user-friendly interface and backlit digital display. Its intuitive design enhances safety and efficiency with features like an automatically locking lid that flashes to signal cycle completion, even in noisy lab environments. Built with durable materials, the HORIZON 12 Flex promises years of maintenance-free operation. Key features include a maximum capacity of 12 tubes, a fully programmable system with three preset cycles for common lab applications, and a unique 2-stage horizontal rotor that ensures optimal sample separation. The centrifuge operates quietly at just 64 dBA and includes a convection-based cooling system to protect samples during heavy workloads. Additional safety features, such as a clear shatter-proof lid and reinforced guard bowl, provide enhanced protection, while the brushless DC motor assures long-term reliability without routine maintenance. Certified to meet US, Canadian, and European safety standards, the HORIZON 12 Flex stands as a robust and efficient choice for any laboratory setting.

[Read More](#)

**SKU:** 00-383-009-000

**Price:** \$2,995.00



## DRUCKER DIAGNOSTICS TRUEBOND VETERINARY CENTRIFUGE FOR BLOOD, URINE, FECAL, AND MICROHEMATOCRIT

The **Drucker Diagnostics TrueBond Veterinary Centrifuge** is an all-in-one solution for veterinary practices, designed to efficiently handle blood, urine, fecal, and microhematocrit samples. This centrifuge comes fully equipped with all necessary accessories, making it a cost-effective choice for veterinary professionals. Included with every TrueBond Centrifuge are six 125 mm tube holders, six universal tube holders suitable for 75-100 mm tubes, six pediatric tube adapters, a microhematocrit reader card, and two inserts for capillary tubes ranging from 40 mm to 75 mm. Its compact design, occupying just 1.35 square feet of counter space, ensures it fits seamlessly into any lab setup. With a maximum speed of 3900 RPM and a capacity of up to six fecal, urine, or blood tubes, or eight microhematocrit tubes, the TrueBond Centrifuge is both powerful and versatile. Users can program up to 10 custom settings via the digital display, tailoring the centrifuge's operation to specific needs. Proudly designed, built, and supported in the USA, this centrifuge also offers the peace of mind of a lifetime warranty on the rotor, underlining its durability and reliability in demanding veterinary environments.

[Read More](#)

**SKU:** 00-376-009-005

**Price:** \$1,045.00



## DRUCKER DIAGNOSTICS HORIZON 6 FLEX FA PROGRAMMABLE ROUTINE CENTRIFUGE

The **Drucker Diagnostics HORIZON 6 Flex FA Programmable Routine Centrifuge** is a versatile and efficient tool designed to meet the demands of routine laboratory work. With its fully digital display and intuitive interface, users can save and recall up to 10 customizable cycle settings, or opt for one of the three expertly crafted presets tailored for chemistry, coagulation, and urine analysis. This centrifuge supports a maximum capacity of six tubes, accommodating sizes up to 16 x 125 mm or 15 mL, and operates with a fixed angle rotor at 45 degrees, ensuring consistent specimen separation. Safety and convenience are paramount, highlighted by the automatic lid lock, shatter-proof clear lid, and reinforced guard bowl. The centrifuge's quiet operation, aided by a brushless DC motor, and its unique convection-based cooling system protect samples under heavy workloads while minimizing maintenance needs. The LED lid lighting system is an innovative feature that prevents forgotten samples and offers clear visual cues for the centrifuge's status.

Certified to meet US, Canadian, and European safety standards, the HORIZON 6 Flex FA combines reliability, efficiency, and user-friendly design, all backed by a two-year warranty and a lifetime guarantee on the rotor, making it an indispensable asset for any laboratory setting.

[Read More](#)

**SKU:** 00-376-009-001

**Price:** \$1,395.00



## DRUCKER DIAGNOSTICS HORIZON 6 FLEX PROGRAMMABLE ROUTINE CENTRIFUGE

The **Drucker Diagnostics HORIZON 6 Flex Programmable Routine Centrifuge** is an exceptional tool designed to meet the varied demands of modern laboratories, particularly where space is limited. With the capability to program and store up to ten custom settings for time, RPM, and g-force, it offers flexibility and efficiency through its user-friendly interface and backlit digital display. The centrifuge features an innovative automatic lid lock with intuitive lid lighting, ensuring safety and ease of use by preventing operation with an open lid and alerting users with flashing lights when cycles are complete. Its robust construction, including a clear, shatter-proof lid and reinforced guard bowl, guarantees sample safety and durability over years of maintenance-free use. The patented 2-stage horizontal rotor design ensures optimal sample separation, while the convection-based cooling system maintains sample integrity even under heavy workloads. With a quiet operation at just 64 dBA, a brushless DC motor, and compliance with international safety standards, the HORIZON 6 Flex is a reliable and advanced solution for routine lab centrifugation needs.

[Read More](#)

**SKU:** 00-376-009-000

**Price:** \$1,995.00



## DRUCKER DIAGNOSTICS HORIZON 24 HIGH-CAPACITY ROUTINE CENTRIFUGE

The **Drucker Diagnostics HORIZON 24 High-Capacity Routine Centrifuge** is a versatile and user-friendly solution for laboratories requiring efficient sample processing. Designed to accommodate a range of tube sizes, including up to 50 mL, this centrifuge offers simplicity through its intuitive set-and-lock controls, ensuring maximum repeatability. Users can easily select from three factory preset cycles or customize settings to fit specific laboratory needs. The centrifuge's lid lights provide clear visual cues regarding cycle status, while the automatic lid lock ensures safety during operation. With features like a patented 2-stage horizontal rotor, brushless DC motor for maintenance-free use, and a convection-based cooling system, the HORIZON 24 maintains optimal sample conditions even under demanding workloads. Certified across multiple safety standards, this centrifuge also boasts a quiet operation at just 59 dBA, making it an ideal choice for busy lab environments. The inclusion of a rotor, tube holders, and a 2-year warranty, along with its affordable pricing, further enhances its appeal as a reliable and robust laboratory instrument.

[Read More](#)

**SKU:** 00-284-009-000

**Price:** \$3,395.00



## DRUCKER DIAGNOSTICS HORIZON 12 COMPACT ROUTINE CENTRIFUGE

The **Drucker Diagnostics HORIZON 12 Compact Routine Centrifuge** is a state-of-the-art device designed to streamline the processing of blood and urine samples, offering exceptional outcomes for patient care. It accommodates up to 16 x 100 mm or 10 mL tubes and features three easily selectable preset cycles based on CLSI guidelines, ensuring optimal performance tailored to specific needs. The centrifuge boasts intuitive LED lid lighting that signals the status of the samples, minimizing errors by preventing forgotten tubes. Its design includes an automatic lid lock for safety and a 2-stage horizontal rotor, ensuring efficient and horizontally separated samples. With a brushless DC motor and a convection-based cooling system, the HORIZON 12 promises years of maintenance-free and cool operation, even under heavy workloads. The centrifuge's quiet operation at just 56 dBA, combined with its clear shatterproof lid for safe observation, makes it a reliable and user-friendly choice for laboratories. Certified to meet US, Canadian, and European safety standards, this centrifuge ensures both durability and compliance, backed by a lifetime rotor warranty.

[Read More](#)

**SKU:** 00-283-009-000

**Price:** \$2,595.00



## **DRUCKER DIAGNOSTICS TRUEBOND SET 3 SET-AND-LOCK VETERINARY CENTRIFUGE FOR BLOOD, URINE, FECAL, AND MICROHEMATOCRIT**

The **Drucker Diagnostics TrueBond Set 3 set-and-lock veterinary centrifuge** is a versatile and compact solution for veterinary laboratories, enabling efficient processing of microhematocrit tubes, as well as blood, urine, and fecal samples. Designed with user-friendly controls, the centrifuge offers preset cycles tailored to common veterinary requirements, ensuring maximum repeatability and ease of use. With a capacity to accommodate up to six blood, urine, and fecal tubes or eight microhematocrits, its fixed-angle rotor and specially designed tube holders cater to various tube sizes, maintaining a 45° angle for optimal results. The TrueBond Set 3 also features intuitive lid lighting to indicate operational status, a convection-based cooling system for sample protection, and an automatic lid lock for safety. Its clear, shatter-proof lid allows for secure sample observation, while the reinforced guard bowl enhances safety and durability. Certified to meet US, Canadian, and European safety standards, the centrifuge is compact and lightweight, fitting easily into tight laboratory spaces and backed by a lifetime warranty on the rotor, making it a reliable and essential tool for veterinary diagnostics.

[Read More](#)

**SKU:** 00-276-009-007

**Price:** \$715.00



## DRUCKER DIAGNOSTICS HORIZON 6 COMPACT ROUTINE CENTRIFUGE

The **Drucker Diagnostics HORIZON 6 Compact Routine Centrifuge** is a versatile and efficient solution for processing blood and urine specimens, particularly in laboratories and point of care facilities where precision and repeatability are paramount. With a capacity to hold up to six tubes and a compact design, it seamlessly fits into any workspace. Users can easily program the centrifuge using its intuitive two-button interface, which allows for the selection of validated cycles or locking into one of three preset cycles recommended by tube manufacturers and CLSI guidelines. This centrifuge is equipped with an innovative LED lid lighting system that indicates the machine's status—off when ready to load, on while running, and flashing at cycle completion—ensuring samples are never forgotten. Safety is a priority, with features like an automatic lid lock preventing operation with an open lid and a reinforced guard bowl for durability. The unique 2-stage horizontal rotor design guarantees optimal sample separation, while the brushless DC motor and convection-based cooling system ensure quiet, maintenance-free operation and sample protection, even under heavy usage. Certified to meet US, Canadian, and European safety standards, the HORIZON 6 promises reliable and safe performance, backed by a lifetime warranty on its rotor.

[Read More](#)

**SKU:** 00-276-009-000

**Price:** \$1,695.00



## DRUCKER DIAGNOSTICS DASH FLEX 24 PROGRAMMABLE STAT CENTRIFUGE

The **Drucker Diagnostics DASH Flex 24 Programmable STAT Centrifuge** is an advanced piece of laboratory equipment designed to optimize efficiency and accuracy in large, high-volume STAT labs. With its innovative features such as LED lid lighting, this centrifuge ensures that samples are never forgotten, thereby reducing turnaround time significantly. The LED lights provide clear visual cues for when the centrifuge is ready, running, or complete, which is especially beneficial in busy and noisy laboratory settings. Its compact design allows it to fit seamlessly in accessioning areas, while its 24-tube capacity facilitates quicker loading and processing. The DASH Flex 24 offers fully programmable settings, enabling users to customize time, RPM, and g-force, or choose from three preset cycles tailored for common STAT requirements. Additionally, its brushless DC motor and unique convection-based cooling system ensure long-lasting, maintenance-free operation and sample integrity under heavy workloads. Safety features like an automatic lid lock, positive lid latching, and a reinforced guard bowl

provide maximum safety and durability, while regulatory approvals guarantee compliance with US, Canadian, and European safety standards. Overall, the DASH Flex 24 is a robust, reliable, and user-friendly centrifuge that enhances laboratory productivity and safety.

[Read More](#)

**SKU:** 00-184-009-000

**Price:** \$4,999.00

## DRUCKER DIAGNOSTICS DASH FLEX 12 PROGRAMMABLE STAT CENTRIFUGE

The **Drucker Diagnostics DASH Flex 12 Programmable STAT Centrifuge** is an advanced laboratory tool designed to optimize efficiency and accuracy in specimen processing. With its innovative features, such as green LED lid lights that clearly indicate the operational status from across the lab, this centrifuge minimizes the risk of samples being overlooked or delayed. Capable of holding up to twelve tubes and accommodating sizes ranging from 75-100 mm, it offers flexibility and convenience in handling various sample volumes. The centrifuge is equipped with a digital display that shows custom parameters and time remaining, and it allows users to program and save up to 10 distinct cycles tailored to specific acceleration, deceleration, speed, g-force, and runtime needs. Its patented 2-stage horizontal rotor ensures optimal sample separation, while the convection-based cooling system protects samples even during intense workloads. Safety is prioritized with features like an automatic lid lock and a clear, shatterproof lid for secure sample observation. The brushless DC motor promises years of maintenance-free operation, and the device adheres to stringent safety regulations, making it a reliable choice for laboratories worldwide. Compact and robust, the DASH Flex 12 is a testament to Drucker Diagnostics' commitment to innovation and quality, offering a lifetime warranty on the rotor for added peace of mind.

[Read More](#)

**SKU:** 00-183-009-000

**Price:** \$4,695.00





## DRUCKER DIAGNOSTICS DASH FLEX 6 DC PROGRAMMABLE STAT CENTRIFUGE

The **Drucker Diagnostics DASH Flex 6 DC Programmable STAT Centrifuge** is a compact yet versatile laboratory tool, perfect for facilities needing efficient and customizable centrifugation without sacrificing bench space. Inspired by the popular DASH Flex 12, this model features a smaller footprint while providing advanced programmable controls and a digital display for easy operation. Users can customize and save up to 10 different cycle parameters, including time, speed, and g-force, tailored to their specific Instrument for Use (IFU) requirements and validation protocols. The centrifuge is pre-configured with the most common STAT settings, capable of spinning up to six tubes with dimensions up to 16 x 100 mm or 10 mL. Its patented swing-out horizontal rotor ensures ideal sample separation and is compatible with tubes from major manufacturers such as Greiner and BD. The DASH Flex 6 prioritizes safety and durability with features like an automatic lid lock, positive lid latching, and a clear shatter-proof lid for safe observation. It also operates with a brushless DC motor for maintenance-free longevity and a convection-based cooling system to protect samples. Certified by US, Canadian, and European safety standards, this centrifuge combines robust performance with user-friendly design, making it a reliable choice for any laboratory setting.

[Read More](#)

**SKU:** 00-176-009-000

**Price:** \$2,695.00



## DRUCKER DIAGNOSTICS DASH APEX 24 HIGH VOLUME STAT CENTRIFUGE

The **Drucker Diagnostics DASH Apex 24 High Volume STAT Centrifuge** is specifically engineered to cater to the demands of large, high-volume STAT laboratories, enhancing efficiency and precision. This advanced centrifuge boasts a compact design, allowing it to fit seamlessly into lab spaces while accommodating up to 24 tubes, facilitating quicker loading and operation. A standout feature is its innovative LED lid lighting system, which provides clear visual cues regarding the centrifuge's status, thereby preventing forgotten samples and reducing turnaround time even in bustling lab environments. The DASH Apex 24 is equipped with an electronic timer offering three easily selectable pre-set cycles, ensuring consistent and repeatable results. Its brushless DC motor guarantees years of maintenance-free operation, while the unique 2-stage horizontal rotor design achieves optimal sample separation. Safety is paramount, with features like an automatic lid lock, positive lid latching, and a clear shatter-proof lid for secure sample observation. Additionally, its convection-based cooling system protects samples during intensive

operations. Certified to meet US, Canadian, and European safety regulations, the DASH Apex 24 is a reliable and efficient choice for labs aiming for high throughput and accuracy.

[Read More](#)

**SKU:** 00-084-009-005

**Price:** \$4,395.00



## DRUCKER DIAGNOSTICS DASH APEX 12 SET-AND-LOCK STAT CENTRIFUGE

The **Drucker Diagnostics DASH Apex 12** Set-and-Lock STAT Centrifuge is an innovative solution designed to enhance efficiency and accuracy in laboratory settings. With its intuitive two-button interface, users can effortlessly select from three pre-set cycles or customize their own validated STAT cycles, significantly minimizing the risk of error and reducing turnaround time. The centrifuge is equipped with LED lid lighting, which provides clear visual cues regarding its operational status, ensuring samples are not forgotten even in the most hectic environments. The DASH Apex 12 can accommodate up to 12 tubes, with sizes up to 16 x 100 mm or 10 mL, and is powered by a brushless DC motor that requires no routine maintenance, guaranteeing years of reliable performance. Its patented 2-stage horizontal rotor ensures optimal separation of samples, while the convection-based cooling system maintains sample integrity under heavy workloads. Safety features such as an automatic lid lock, positive lid latching, a clear shatterproof lid, and a reinforced guard bowl are integrated to protect both users and samples. Certified to meet US, Canadian, and European safety standards, this centrifuge is both a robust and economical choice, supplied with all necessary components, including a lifetime warranty on the rotor, ensuring peace of mind for laboratory professionals.

[Read More](#)

**SKU:** 00-083-009-003

**Price:** \$3,395.00



## DRUCKER DIAGNOSTICS 642M 12V DC MOBILE CENTRIFUGE

The **Drucker Diagnostics 642M 12V DC Mobile Centrifuge** is an innovative solution tailored for home health nurses who require efficient and reliable blood sample processing on the go. This compact centrifuge, powered by a vehicle's 12V electrical system, is engineered to provide optimal pre-analytic horizontal-sample-processing, aligning with the recommendations of all blood collection tube manufacturers. Utilizing patented technologies, the 642M accelerates blood separation by 33%, significantly reducing wait times and allowing healthcare professionals to maximize their time between patient appointments. The centrifuge boasts a variety of features, including a preset electronic timer for ideal gel separations, a capacity for up to six tubes, a maintenance-free 12 Volt DC motor, and a double-encased design for quiet operation with minimal vibration. Its unique convection-based cooling system ensures sample integrity even under heavy workloads, while the compact footprint requires only about one square foot of counter space. Additional safety and convenience features include an automatic lid lock, a clear shatter-proof lid for safe observation, and a reinforced guard bowl for enhanced durability. With a maximum speed of 3400 RPM, rapid acceleration, and deceleration times, the 642M is a reliable and efficient tool for any mobile health professional, backed by a lifetime warranty on its rotor.

[Read More](#)

**SKU:** 00-077-009-000

**Price:** \$745.00



## DRUCKER DIAGNOSTICS DASH APEX 6 COMPACT STAT CENTRIFUGE

The **Drucker Diagnostics DASH Apex 6 Compact STAT Centrifuge** is a game-changer for laboratories striving to enhance efficiency and accuracy in healthcare delivery. Designed to meet the demands of fast-paced lab environments, the Apex 6 effectively halves STAT turnaround times, providing high-quality plasma samples in just minutes. Its compact design allows for seamless integration into any lab setup, fitting easily onto tabletop areas for immediate sample processing. The centrifuge's user-friendly two-button interface eliminates the need for complex programming, enabling quick and hassle-free operation straight out of the box. Capable of accommodating up to six tubes, each up to 16 x 100 mm or 10 mL, the Apex 6 does away with the need for costly unique holders. Its advanced features, such as a patented two-stage horizontal rotor and a convection-based cooling system, ensure optimal sample integrity and safety. Certified to meet safety regulations across the US, Canada, and Europe, the Apex 6 also boasts a brushless DC motor for maintenance-free

operation and an automatic lid lock for enhanced security. With preset cycles that cater to various processing times and speeds, the DASH Apex 6 promises to improve workflow and deliver reliable results, making it an indispensable tool for any lab at the forefront of healthcare.

[Read More](#)

**SKU:** 00-076-009-004

**Price:** \$3,195.00



## DRUCKER DIAGNOSTICS TRUEBOND DUET COMPACT DUAL ROTOR VETERINARY CENTRIFUGE

The **Drucker Diagnostics TrueBond Duet Compact Dual Rotor Veterinary Centrifuge** is an essential piece of equipment for veterinary laboratories, providing a versatile and compact solution for processing various samples. This innovative centrifuge is equipped with two removable rotors, allowing seamless transitions between spinning 1.5 - 2 mL blood or urine tubes and accommodating up to twelve 40 mm microhematocrit tubes. Its user-friendly design includes five easily selectable preset cycles, ensuring consistent and reliable results with every use, making it highly valuable for any veterinary practice. The TrueBond Duet's unique convection-based cooling system ensures cool operation, protecting samples even during heavy workloads. Safety features include an automatic lid lock, a clear shatter-proof lid for safe observation, and a reinforced guard bowl, all of which contribute to its reliability and durability. Certified to meet US, Canadian, and European safety standards, this centrifuge is compact in size and lightweight, making it a practical choice for any laboratory setting. With a lifetime warranty on the rotor, the TrueBond Duet promises long-term performance and peace of mind for veterinary professionals.

[Read More](#)

**SKU:** 00-056-009-000

**Price:** \$1,895.00



## DRUCKER DIAGNOSTICS DASH APEX 4 COMPACT STAT CENTRIFUGE

The **Drucker Diagnostics DASH Apex 4 Compact STAT Centrifuge** is an innovative solution for laboratories seeking efficiency and convenience without compromising on performance. Weighing just 4.3 lbs and with dimensions of 7 x 7 x 7 inches, it is the smallest model within the DASH STAT centrifuge family, making it perfect for labs with limited bench space. Its user-friendly two-button interface allows for easy selection of up to three preset cycles or customization for consistent results. The DASH Apex 4 can accommodate up to four tubes, with a maximum capacity of 5 mL, and features a brushless DC motor that requires no routine maintenance, ensuring longevity and reliability. The centrifuge's fixed angle rotor is compatible with all common tube manufacturers, while its unique convection-based cooling system maintains sample integrity during intensive use. Safety is paramount, with features like an automatic lid lock, positive lid latching, and a clear shatter-proof lid for secure operation and observation. Certified to meet US, Canadian, and European safety standards, the DASH Apex 4 is a compact powerhouse that delivers precision and safety in a streamlined design.

[Read More](#)

**SKU:** 00-053-009-000

**Price:** \$2,995.00



## DRUCKER DIAGNOSTICS 642E SINGLE- SPEED ELECTRONIC CENTRIFUGE

The **Drucker Diagnostics 642E Single-Speed Electronic Centrifuge** is a meticulously engineered device that stands out for its optimal pre-analytic sample processing capabilities, making it an indispensable tool in medical laboratories and doctor's offices. With its compact design, utilizing just about one square foot of counter space, the 642E is perfect for environments where space is at a premium. Its patented rotor technology significantly reduces blood separation processing times by 33%, thereby enhancing efficiency without compromising sample integrity. The centrifuge features a brushless AC motor, ensuring years of maintenance-free operation, and operates quietly due to its double-encased motor and stability mount design. Safety and ease of use are prioritized with features like an automatic lid lock, a safety lid switch, and a clear shatter-proof lid for safe observation. The centrifuge is capable of handling up to six tubes, with a maximum speed of 3380 RPM, and its convection-based cooling system keeps samples cool even under heavy workloads. Certified to US, Canadian, and European safety standards, the 642E's robust build, including a reinforced guard bowl and positive lid latching, ensures maximum safety and durability. Its lifetime warranty on

the rotor further underscores its reliability and appeal as the world's best-selling clinical laboratory centrifuge.

[Read More](#)

**SKU:** 00-018-009-001

**Price:** \$755.00



## DRUCKER DIAGNOSTICS 642B SINGLE-SPEED CENTRIFUGE

The **Drucker Diagnostics 642B Single-Speed Centrifuge** is a cost-effective solution for laboratories seeking reliable and efficient sample processing. This model, renowned for its affordability, operates at a consistent 1600 xg, utilizing Drucker's patented horizontal rotor technology to deliver superior sample separations. Its advanced cooling system ensures sample integrity even under heavy workloads, while the brushless AC motor promises years of maintenance-free operation. The 642B accommodates up to six tubes of varying sizes, offering versatility for a wide range of applications with a mechanical timer that ranges from 1 to 30 minutes. Its compact design, occupying just one square foot of counter space, is complemented by features like a clear shatterproof lid for safe observation, reinforced side walls for added safety, and a quiet operation due to its double-encased motor. Certified to meet US and Canadian safety standards, the centrifuge also includes a positive lid latching system to secure samples during operation. With a lifetime warranty on the rotor, the 642B ensures peace of mind and dependable performance for any laboratory setting.

[Read More](#)

**SKU:** 00-016-009-001

**Price:** \$735.00



## DRUCKER DIAGNOSTICS 614B FIXED ANGLE SINGLE-SPEED CENTRIFUGE

The **Drucker Diagnostics 614B Fixed Angle Single-Speed Centrifuge** is an efficient and cost-effective solution for medical settings that require reliable specimen processing. It offers a single g-force setting of 1100 xg and a mechanical timer adjustable from 1 to 30 minutes, simplifying operations by eliminating the need for complex settings. This centrifuge ensures consistent and error-free chemistry and coagulation spins, making it ideal for use in doctor's offices, small clinics, and laboratories where on-site analysis is performed. However, while its fixed-angle rotor is effective for up to six tubes, it's important to note that gel layers formed may be vulnerable to rupture during transport, necessitating re-centrifugation, which can negatively affect analytical results. For specimens that will be transported, horizontal centrifugation models like the Model 642E or HORIZON series are recommended. The 614B's compact design, occupying just one square foot of counter space, features a brushless AC motor for long-lasting, maintenance-free operation, quiet and cool performance, and multiple safety features, including a positive lid latching system and reinforced side walls. With regulatory approvals and a lifetime warranty on the rotor, the 614B provides a reliable and secure solution for various laboratory applications.

[Read More](#)

**SKU:** 00-004-009-001

**Price:** \$645.00



## STATSPIN EXPRESS 4 CENTRIFUGE

The **StatSpin Express 4 high-speed horizontal centrifuge** revolutionizes laboratory efficiency by providing rapid turnaround times and enhancing productivity, particularly in Lean laboratory settings. Positioned adjacent to any STAT analyzer, it ensures that serum or plasma gel tubes are ready in mere minutes, thanks to its horizontal spin that creates a flat gel barrier. The device is user-friendly, featuring only three settings—3, 5, or 10 minutes—eliminating the need for speed adjustments or complicated settings. Safety is paramount, with a lid that remains locked until the rotor stops completely and an alarm that signals the end of the cycle. This allows for immediate processing of samples as they arrive, removing the inefficiencies of batching and waiting. Capable of handling up to eight samples simultaneously, the Express 4 significantly boosts throughput. In just three minutes, it delivers high-quality serum or plasma separation, making it ideal for rapid separation of gel tubes, creating flat gel barriers, preparing samples for STAT chemistry, and facilitating single-piece flow for swift sample processing. Key features include its compact, fast, and powerful design for improved productivity,

ease of use, universal power supply for global compatibility, maintenance-free brushless motor, and patented rotor design.

[Read More](#)

**SKU:** AM-SS-EX4

**Price:** \$2,890.00



## SEROFUGE 2002 CENTRIFUGE

The **Clay Adams Sero-Fuge, Model 2002 Centrifuge**, is a versatile and highly efficient system designed to meet the rigorous demands of blood grouping, cross matching, typing, and other cell washing procedures. This advanced centrifuge from Clay Adams boasts a 12-place head rotor and operates at two speeds, making it adaptable for various laboratory needs. Its superior features include a locking lid for enhanced safety, digital controls and displays for precise operation, and a gasket seal to ensure optimal performance. The dual-speed operation allows for settings of either 2400rpm/500xG or 3400rpm/1049xG, providing flexibility depending on the specific requirements of the procedure. Additionally, the Sero-Fuge 2002 includes an unbreakable 12-place plastic rotor designed to accommodate blood tubes up to 12.7 mm in diameter and 75 mm in length, ensuring durability and reliability. This centrifuge stands out as a top choice for laboratories seeking high performance and versatility in their cell washing and blood testing processes.

**LEAD TIME: 6 Weeks**

[Read More](#)

**SKU:** AM-SEROFUGE-2002

**Price:** \$2,995.00

# INDEX

## D

- Drucker Diagnostics 614B Fixed Angle Single-Speed Centrifuge 17
- Drucker Diagnostics 642B Single-Speed Centrifuge 16
- Drucker Diagnostics 642E Single-Speed Electronic Centrifuge 15
- Drucker Diagnostics 642M 12V DC Mobile Centrifuge 13
- Drucker Diagnostics DASH Apex 12 Set-and-Lock STAT Centrifuge 12
- Drucker Diagnostics DASH Apex 24 High Volume STAT Centrifuge 11
- Drucker Diagnostics DASH Apex 4 Compact STAT Centrifuge 15
- Drucker Diagnostics DASH Apex 6 Compact STAT Centrifuge 13
- Drucker Diagnostics DASH Coag 2-3-5 Minute STAT Coag Centrifuge 2
- Drucker Diagnostics DASH Flex 12 Programmable STAT Centrifuge 10
- Drucker Diagnostics DASH Flex 24 Programmable STAT Centrifuge 9
- Drucker Diagnostics DASH Flex 6 DC Programmable STAT Centrifuge 11
- Drucker Diagnostics HORIZON 12 Compact Routine Centrifuge 7
- Drucker Diagnostics HORIZON 12 Flex Programmable Routine Centrifuge 4
- Drucker Diagnostics HORIZON 24 Flex Programmable Routine Centrifuge 3
- Drucker Diagnostics HORIZON 24 High-Capacity Routine Centrifuge 7
- Drucker Diagnostics HORIZON 6 Compact Routine Centrifuge 9
- Drucker Diagnostics HORIZON 6 Flex FA Programmable Routine Centrifuge 5
- Drucker Diagnostics HORIZON 6 Flex Programmable Routine Centrifuge 6
- Drucker Diagnostics SERO 12 Blood Banking Centrifuge 3
- Drucker Diagnostics TrueBond Duet Compact Dual Rotor Veterinary Centrifuge 14
- Drucker Diagnostics TrueBond Set 3 Set-and-lock veterinary centrifuge for blood, urine, fecal, and microhematocrit 8
- Drucker Diagnostics TrueBond Veterinary centrifuge for blood, urine, fecal, and microhematocrit 5

## S

- Serofuge 2002 Centrifuge 18
- StatSpin Express 4 Centrifuge 17