

TABLE OF CONTENTS

Communications 2

 C-Com Lite Headset System 3



COMMUNICATIONS

Effective medical communication is vital for enhancing collaboration among healthcare professionals, improving patient outcomes, and ensuring safety in clinical environments. Clear communication fosters teamwork, particularly in high-pressure situations such as surgeries or emergency care, where miscommunication can lead to critical errors. By utilizing advanced communication tools, medical teams can streamline information sharing, coordinate patient care more efficiently, and reduce the likelihood of misunderstandings, ultimately leading to better patient experiences and satisfaction.

In our Medical Communications category, we offer innovative solutions designed to facilitate seamless communication within clinical teams. One standout product is the C-Com Lite Headset System, which enhances audio clarity and allows for hands-free communication in environments like Cath Labs and Operating Rooms. This system is specifically engineered to improve collaboration among team members, ensuring that critical instructions and information are conveyed clearly and promptly. By incorporating such advanced communication tools, healthcare facilities can significantly bolster their operational efficiency and enhance the overall quality of care provided to patients.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

[REFURBISHED](#)

C-COM LITE HEADSET SYSTEM



The **C-Com Lite Headset System** enhances communication within clinical teams, particularly in Cath, EP, IR, and Hybrid OR environments. With whisper-sensitive individual microphones, the system ensures flawless audio fidelity, enabling precise instructions and collaborative communication to be effortlessly shared among all team members. This not only keeps everyone engaged but also fosters a calmer lab environment, which directly benefits patient care and outcomes. The C-Com Lite Headset System boasts impressive specifications, including a 16-hour battery life, the capability for four simultaneous talkers and listen pathways, secure transmission, and a noise-canceling microphone. It's designed for easy integration with both new and existing sound or intercom systems, making it a versatile addition to any clinical setting. The system dimensions are 9.5" x 5" x 3" and it weighs just 1lb, ensuring portability and convenience. Operating on a 2.4-2.4835 GHz FHSS frequency, it offers up to a 300 ft operating range. Additionally, it comes with a charging station featuring four active and four reserve charging bays, ensuring the headset is always ready to go. The C-Com Lite Headset System is powered by a removable, rechargeable lithium-ion battery with a low battery indicator, making it a reliable choice for continuous clinical use.

[Read More](#)

SKU: AM-C-COM-LITE

INDEX

C

C-Com Lite Headset System 2