

TABLE OF CONTENTS

ECG	2
Trackmaster TMX 425 Treadmill	3
Trackmaster TMX 428 Treadmill	3
GE MAC 5500 HD ECG	4
SunTech Tango M2 Blood Pressure (NEW) [To be paired with Norav Stress]	4
Norav Stress System Computer Tower	5
Norav Stress System, Wireless, Win 10, with Treadmill	5
Mortara ELI 280 ECG Machine	5
GE Mac 2000 ECG Machine, Basic 12-lead Resting ECG System	6
Mortara ELI 250 ECG Machine	7
Burdick 8500 ECG Machine	7
Welch Allyn Schiller AT-2 Plus	7
Edan SE-1515 PC-based ECG DX12 w/ ECG Receiving box	8
Welch Allyn CP150	8
Atria 6100	9
GE Mac 5000 EKG	9
Burdick 8300	10
Edan SE-601 Series ECG	10
GE MAC 1200	11
Edan SE-1200 and SE-1200 Express ECG	11



ECG

Cutting-edge electrocardiogram (ECG) technology helps hospitals and medical practices maintain optimum heart health for their patients. Auxo Medical provides a wide scope of popular ECG machines, including [Edan](#), [GE](#), and [Burdick](#).

These easy-to-use ECG machines allow your medical facility to operate smoothly and efficiently. Furthermore, our machines offer real time ECG waveform displays, compact and lightweight designs, one-touch operations, and alphanumeric keyboards. In addition, many of Auxo Medical's ECG machines offer powerful data analysis functions to promote accuracy and maximum decision-making efficiency for any medical practitioner.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

REFURBISHED



TRACKMASTER TMX 425 TREADMILL

The **Trackmaster TMX 425 Treadmill** is a robust, high-capacity stress-testing treadmill designed for precision and durability in clinical and research settings. Equipped with a heavy-duty 2.2 HP AC inverter drive and a double-sided, lubricated, cushioned running deck measuring 22" x 63", it provides a smooth, low-impact walking surface at a low profile of 7.5 inches from the floor. The treadmill supports a user capacity of up to 500 lbs and features adjustable elevation from 0-25%, with speeds ranging from 0.5 to 12 MPH on 220V power or 0.5 to 10 MPH on 110V. Constructed with 10-gauge steel side rails, a 12-gauge steel motor pan, and epoxy powder coating in jet-black with a white finish, it combines strength with stylish durability. Standard safety side rails, an RS232 interface, and compatibility with the QRS system make it a reliable choice for rigorous testing environments. The treadmill requires a level 33" x 88" surface for optimal use and offers advanced features including auto speed reset, ensuring precision and safety during high-stress evaluations.

[Read More](#)

SKU: AM-TMX-425

Price: Request Quote



TRACKMASTER TMX 428 TREADMILL

The **Trackmaster TMX 428 Treadmill** is a premier stress-testing system engineered for precision, safety, and durability. Featuring the widest industry-standard 22-inch running belt, it provides users with unmatched confidence during high-intensity tests, supported by the patented MasterTrack v-guide that ensures perfect belt alignment and eliminates lateral movement. Mounted on shock-absorbing isolators, the deck offers enhanced impact attenuation for patient comfort and safety. Its smooth, quiet operation is powered by a high-performance 4 HP continuous PWM DC brushless servo motor, capable of reaching speeds up to 15 MPH with self-calibrating speed and elevation controls. The low-profile, cushioned, reversible deck, adjustable elevation from 0-25%, and emergency stop options—including a prominent right or left hand button and tether—make this treadmill both reliable and safe. Built with an all-steel chassis, minimal maintenance requirements, and advanced connectivity ports, the TMX 428 exemplifies cutting-edge engineering designed for rigorous stress testing in clinical and research environments.

[Read More](#)

SKU: AM-TMX-428

Price: Request Quote



GE MAC 5500 HD ECG

The **GE MAC 5500 HD Resting ECG System** is a high-end ECG system designed for speed, accuracy, and comprehensive disease management. Built upon GE's legacy of ECG innovation, the MAC 5500 HD offers advanced algorithms and networking capabilities, making it a powerful tool for clinicians. Key features include improved triaging for Acute Coronary Syndrome (ACS) detection, configurable critical values to accelerate treatment, and enhanced pacemaker detection for accurate interpretation of paced rhythms. The system also streamlines workflow with barcode capabilities for order and ADT downloads, aiming to boost efficiency across healthcare enterprises. Ultimately, the MAC 5500 HD facilitates connected ECG workflows, helping to save time, reduce misdiagnoses, and minimize errors in patient data and billing, allowing clinicians to focus on delivering optimal patient care.

[Read More](#)

SKU: AM-MAC5500

Price: Request Quote



SunTech Tango M2

Blood Pressure (NEW)

To be paired with Norav Stress

SUNTECH TANGO M2 BLOOD PRESSURE (NEW) [TO BE PAIRED WITH NORAV STRESS]

Includes: 2 Cuffs w/mic, hose, power cable

[Read More](#)

SKU: AM-ST-TANGO-M2

Price: Request Quote



Norav Stress System Computer Tower

NORAV STRESS SYSTEM COMPUTER TOWER

The Norav Stress System Computer Tower package is a comprehensive solution perfect for Norav Stress testing. This bundle includes a new Windows 10 computer tower for maximum compatibility. Accompanying the tower are essential peripherals: printer, monitor, keyboard, and mouse. Additionally, the package contains a sturdy cart for easy mobility. All systems come preloaded with AIM software. However, it's important to note that **setup is not included**.

[Read More](#)

SKU: AM-NORAV-CT

Price: Request Quote



Norav Stress System, Wireless, Win 10, with Treadmill

NORAV STRESS SYSTEM, WIRELESS, WIN 10, WITH TREADMILL

The Norav Stress System represents a state-of-the-art solution for cardiovascular assessment, seamlessly integrating with Windows 10 through a wireless module and including a robust treadmill designed for clinical settings. This system is not only user-friendly but also equipped with advanced diagnostic tools that connect easily to any PC via USB. It boasts a comprehensive suite of features aimed at enhancing diagnostic accuracy and patient safety, including a standard 12-lead detachable cable, built-in defibrillation protection, and automatic arrhythmia detection. The free software updates for the lifetime of the product ensure that practitioners always have access to the latest advancements. Additionally, the system supports a variety of protocols and offers automatic control of treadmills and optional ergometers. With capabilities like Heart Rate Recovery Score, dynamic ST monitoring, and detailed PDF report exports, the Norav Stress System is a cost-effective yet powerful option for both office and hospital use. **Please note, the setup does not include a Computer Tower, and installation is not included.**

[Read More](#)

SKU: AM-NORAV-T

Price: Request Quote



MORTARA ELI 280 ECG MACHINE

The **Mortara ELI™ 280 electrocardiograph** represents a pinnacle of portability and high-level functionality, designed to thrive in the demanding conditions of high-volume ECG environments. This compact, lightweight device features a 10.1" touchscreen high-resolution color LCD, enhancing the display of up to 12 leads of ECG data for efficient review and analysis. Its user-friendly interface and widescreen layout simplify navigation and patient data entry, while ensuring clear visibility of electrode connections. The ELI 280 meets the latest FAA standards for EKG report submissions, making it ideal for aviation medical examinations. Among its standout features is the "Best 10" function which selects the optimal 10 seconds of ECG data from memory, thus reducing the need for repeat tests due to artifacts like patient movement. Equipped with a Wired Module for enhanced connectivity, the ELI 280 supports PDF outputs and seamless integration with EHR/MIS systems through DICOM®. Its advanced data capture capability at 1,000 sps per channel, and exceptional pacemaker spike detection at 40,000 sps per channel, set new standards in ECG precision and reliability. This robust device comes complete with essential accessories such as a 10 lead patient cable, electrodes, and paper, ensuring it is ready for immediate deployment in diverse medical settings.

[Read More](#)

SKU: AM-ELI280

Price: Request Quote



GE MAC 2000 ECG MACHINE, BASIC 12-LEAD RESTING ECG SYSTEM

The **GE Mac 2000 ECG** system provides proven diagnostic support enabling clinicians to make fast, accurate diagnoses. Its flexible connectivity solutions offer clinicians the data they need to help improve patient outcomes and streamline workflow throughout the continuum of cardiac care.

The GE Mac 2000 Workflow enhancements with connectivity to EMRs, CS Diagnostic System and the MUSE* ECG Management systems, or send reports to a shared directory on a PC

GE Mac 2000 - Busy clinical practices need streamlined data connectivity, without the time-consuming complexity. That's why the GE Mac 2000 is designed to connect with your existing clinical systems, keeping you in touch with your data, your patients, and your practice.

[Read More](#)

SKU: AM-MAC2000-B

Price: Request Quote

MORTARA ELI 250 ECG MACHINE

Compact and lightweight, the ELI 250c Electrocardiograph provides comprehensive functionality in an easy-to-use, portable device. Multiple connectivity options ensure patient data transfers easily from the ECG to the location of your choice. Experience the benefits of bidirectional communication via internal modem, LAN or wireless LAN with Mortara's E-Scribe™ and Athena products, as well as third-party EMR, HIS, and PACS systems via XML, PDF and DICOM®. Best 10 technology helps you capture the best 10 seconds of ECG data from system memory, reducing the need for repeat ECGs. Wireless data acquisition helps increase workspace flexibility by gathering 12-lead data without the need for a tethered cable.


[Read More](#)**SKU:** AM-ELI250**Price:** Request Quote

BURDICK 8500 ECG MACHINE

The Burdick 8500 ECG / EKG shows you what you need, comfortably. Burdick 8500 ECG equipment has a friendly user interface, multi-angle 7-inch widescreen, and color display that makes reading waveform data easy, and protocol choices streamline your workflow.

[Read More](#)**SKU:** AM-B8500**Price:** Request Quote

WELCH ALLYN SCHILLER AT-2 PLUS



SCHILLER's AT-2 light offers the economics of a low cost unit with many of the high end features needed to care for your patients. The AT-2 light is a 12-channel ECG system offering battery operation of up to 6 hours, full size (8 1/2 x 11) reports to fit patient records, LCD screen and keyboard for easy entry of patient demographics and interpretation.

[Read More](#)

SKU: AM-WA-SCH-AT2+

Price: Request Quote



EDAN SE-1515 PC-BASED ECG DX12 W/ ECG RECEIVING BOX

The SE-1515 integrates high performance PC software with a compact yet powerful EKG machine. By combining reliable algorithms with rich outputs, this all-in-one workstation can help maximize your daily workflow.

[Read More](#)

SKU: AM-SE-1515-DX12

Price: Request Quote

WELCH ALLYN CP150



The Welch Allyn CP 150™ Resting Electrocardiograph (ECG) with optional spirometry combines the best traditional box ECG features like full-size printouts and clinical decision with a simple touchscreen interface, fast one-button operation, and flexible connectivity options to help you keep your staff productive and your tests accurate. And with the optional spirometry module, the CP 150 even lets you perform two reimbursable tests with just one device. The end result is a device that will help your facility save time, improve workflow, and efficiently manage patient information.


[Download Brochure](#)

[Read More](#)

SKU: AM-WA-CP150

Price: Request Quote

ATRIA 6100



The Atria 6100 is a practical, portable ECG / EKG machine for hospitals and cardiology clinics. The Atria 6100 portable ECG helps streamline workflow with a friendly user interface and clear, readable color display for easy reading of ECG waveform data.

[Download Brochure](#)

[Read More](#)

SKU: AM-ATRIA6100

Price: Request Quote



GE MAC 5000 EKG

MAC® 5000, GE's premier ECG system, delivers advanced disease management capabilities through its vast collection of proven algorithms. The MAC 5000 system offers the sophistication required for advanced ECG applications, while its ease of use extends this level of performance to the broadest range of users possible. And it's part of the complete GE suite of networked, non-invasive testing solutions designed to maximize patient throughput and department productivity. Add the stress option, and you combine the industry's leading high performance resting ECG system with comprehensive, equally advanced exercise-testing capabilities in a single, compact system.

[Download Brochure](#)

[Read More](#)

SKU: AM-MAC5000

Price: Request Quote

BURDICK 8300

The Burdick 8300 is a practical, no-nonsense ECG device. The Burdick 8300 ECG / EKG has a friendly user interface and clear, readable 5.7-inch monochrome display to make reading data easy, with protocol choices designed to streamline your work flow.

[Download Brochure](#)

[Read More](#)

SKU: AM-B8300

Price: Request Quote



EDAN SE-601 SERIES ECG



Six-channel ECG SE-601 series : SE-601A, SE-601B, SE-601C, cost-effective solutions including hospital wards, routine check, home visit and emergency cases.

[Download Brochure](#)

[Read More](#)

SKU: N/A

Price: Request Quote



GE MAC 1200

The MAC* 1200 digital, 12-lead electrocardiograph system offers comprehensive ECG solutions with practical features that help meet the needs of hospitals, clinics, office based practices and clinical trials. Its advanced algorithm capabilities, seamless connectivity to the CS and MUSE* Cardiology Information Systems and easy-to-use features provide the quality and convenience users require.

[Download Brochure](#)

[Read More](#)

SKU: AM-MAC1200

Price: Request Quote



EDAN SE-1200 AND SE-1200 EXPRESS ECG

Excellent cost performance with built-in high-resolution thermal printer. You can record 12 ECG traces on A4 size and US Letter-size paper.

[Edan SE-1200 Brochure](#)

[Edan SE-1200 Express Brochure](#)

[Read More](#)

SKU: AM-SE-1200EX

Price: Request Quote

INDEX

A

Atria 6100 8

B

Burdick 8300 9

Burdick 8500 ECG Machine 6

E

Edan SE-1200 and SE-1200 Express ECG 11

Edan SE-1515 PC-based ECG DX12 w/ ECG Receiving box 7

Edan SE-601 Series ECG 10

G

GE MAC 1200 10

GE Mac 2000 ECG Machine, Basic 12-lead Resting ECG System 5

GE Mac 5000 EKG 9

GE MAC 5500 HD ECG 3

M

Mortara ELI 250 ECG Machine 6

Mortara ELI 280 ECG Machine 5

N

Norav Stress System Computer Tower 4

Norav Stress System, Wireless, Win 10, with Treadmill 4

S

SunTech Tango M2 Blood Pressure (NEW) [To be paired with Norav Stress] 3

T

Trackmaster TMX 425 Treadmill 2

Trackmaster TMX 428 Treadmill 2

W

Welch Allyn CP150 8

Welch Allyn Schiller AT-2 Plus 7