

TABLE OF CONTENTS

ENT	2
Mindray V12 Monitor Only	3
Mindray V12 V21 Nellcor Superpod Module	3
Alaris 8120 PCA Syringe Pump	4
Midmark 498 ENT Cabinet With Vacuum	6
Mobile ENT Cart	6



ENT

Auxo Medical provides the necessary cabinetry for ears, nose, and throat personnel. Specifically designed for easy accessibility and portability, these cabinets will help keep your instruments and equipment sterile stay organized.

Auxo Medical offers cabinets with lockable, mobile casters that have the capacity to move and support all necessary instruments. Built from solid steel, with features including an endosoak system, fiber optic lights, suction and pressure pumps, these durable, state-of-the-art ENT carts. It offers everything your facility needs to maintain operations and top quality patient care. Auxo Medical understands how important mobility is when dealing with ENT services and procedures.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

RECONDITIONED

REFURBISHED

MINDRAY V12 MONITOR ONLY

The **Mindray V12 Monitor** is a versatile, modular patient monitoring platform featuring a 12.1-inch color TFT touchscreen display with up to 8 waveforms, designed to adapt to diverse clinical needs. Its flexible architecture allows for easy addition or removal of modules to customize monitored parameters based on patient acuity or monitoring requirements, making it suitable for a wide range of hospital settings. The V12 supports integration with Panorama Central Stations and Nurse Call Systems, and offers advanced features such as IV drug and hemodynamic calculations, a modern user interface, and easy-to-maintain interchangeable system elements. With a compact design—measuring approximately 13.7 inches high, 9.4 inches wide, and 7.5 inches deep—and weighing just over 8.9 pounds (without batteries), it provides reliable, high-performance monitoring in a lightweight, space-efficient package. Please note that a Mindray Parameters Module must be purchased separately for patient monitoring.

[Read More](#)

SKU: AM-MR-V12

Price: \$2,195.00





MINDRAY V12 V21 NELLCOR SUPERPOD MODULE

The **Mindray V12 V21 Nellcor Superpod Module** is a refurbished, highly versatile monitoring accessory designed to enhance the capabilities of Mindray V Series monitors. This module provides comprehensive patient monitoring with ECG and respiration assessment, temperature measurement, dual invasive blood pressure (IBP), Nellcor SpO₂, and non-invasive blood pressure (NIBP). Compatible with V12 and V21 models, it offers a flexible, integrated solution for adult, pediatric, and neonatal care. Please note that this listing includes only the module; accessories are available separately for complete setup. OEM Part Number: 0998-00-1802-0102A.

[Read More](#)

SKU: 0998-00-1802-0102A

Price: \$795.00



ALARIS 8120 PCA SYRINGE PUMP

The **CareFusion Alaris 8120 PCA Syringe Pump** is an advanced infusion device designed for safe and effective patient-controlled analgesia (PCA) management. Integrated within the Alaris® System, it combines programming safety features, barcoding capabilities, and continuous respiratory monitoring into a comprehensive pain control solution. The PCA module offers quick access to total dose history and trending data, enabling healthcare providers to analyze patient responses and optimize pain therapy. Notably, it can automatically pause PCA infusions and deactivate dose requests if the patient's oxygen saturation (SpO₂) or exhaled CO₂ (EtCO₂) levels fall below hospital-defined thresholds when used with the PCA module. The pump allows for different dose ranges for the same drug based on patient diagnosis (e.g., acute vs. chronic pain) and supports barcoding to ensure the accuracy of drug, patient, and clinician identification. Its intuitive interface simplifies programming, reducing training time, while the Guardrails® Suite MX software enables customizable alarm limits and trend monitoring. When combined with the Auto-ID module, it provides reliable barcoding for the right drug, patient, and clinician, enhancing safety and workflow efficiency in pain management.

Warranty: 6 Months

[Read More](#)

SKU: AM-Alaris-8120

Price: \$2,261.00



MIDMARK 498 ENT CABINET WITH VACUUM

The Midmark 498 power procedure center provides excellent storage space, holds instruments, provides strong aspiration (suction), and adjustable pressure. A must-have piece of furniture for the medical professional. The major components of the unit consist of a drawer assembly, power on/standby switch, vacuum on/off switch, pressure switch, regulator, pressure gauge, front receptacle, rear receptacle, vacuum/pressure foot switch, pressure pump assembly, and vacuum pump assembly.

[Read More](#)

SKU: AM-M498





MOBILE ENT CART

The Mobile ENT Cart is an all-steel, lockable office on heavy-duty casters. Designed to hold all instruments and equipment needed for ENT, diagnostics, and minor procedures. Solid-steel construction insures the cart is ready when you need it.

[Read More](#)

SKU: AM-03-7570

INDEX

A

Alaris 8120 PCA Syringe Pump 3

M

Midmark 498 ENT Cabinet With Vacuum 5

Mindray V12 Monitor Only 2

Mindray V12 V21 Nellcor Superpod Module 3

Mobile ENT Cart 6